

PCTWORLD INTELLECTUAL PROPERTY ORGANIZATION
International Bureau

INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification 5 : C12Q 1/68	A2	(11) International Publication Number: WO 93/03187 (43) International Publication Date: 18 February 1993 (18.02.93)
(21) International Application Number: PCT/US92/06617 (22) International Filing Date: 28 July 1992 (28.07.92) (30) Priority data: 738,800 31 July 1991 (31.07.91) US (71) Applicant: AMOCO CORPORATION [US/US]; Mail Code 1907, Patents and Licensing Department, P.O. Box 87703, Chicago, IL 60680-0703 (US). (72) Inventors: PARODOS, Kyriaki ; 30 Royal Crest Drive, Apartment 9, Marlborough, MA 01752 (US). McCARTY, Janice, Marie ; 10 Donna Terrace, Hyde Park, MA 02136 (US).		(74) Agents: LADD, Robert, G. et al.; Amoco Corporation, Patents and Licensing Department, Mail Code 1907, P.O. Box 87703, Chicago, IL 60680-0703 (US). (81) Designated States: JP, European patent (AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LU, MC, NL, SE). Published <i>Without international search report and to be republished upon receipt of that report.</i>
(54) Title: NUCLEIC ACID PROBES FOR THE DETECTION OF SHIGELLA (57) Abstract The invention relates to methods of detection of bacteria of the genus <i>Shigella</i> and/or Enteroinvasive <i>E. coli</i> (EIEC) by use of a set of nucleic acid probes. The invention further relates to a set of <i>Shigella</i> specific chromosomal sequences and fragments and to probes derived from the <i>Shigella</i> specific fragments. Additionally, probes were derived from a sequence from the <i>Shigella</i> <i>ompA</i> gene. In particular, a series of probes, each approximately 40 nucleotides in length, were designed having specificity for <i>Shigella</i> or for <i>Shigella</i> and Enteroinvasive <i>E. coli</i> , and having utility in nonisotopic test formats which require amplification to achieve high sensitivity. Specific hybridization probe sets which are capable of detecting substantially all clinically significant serotypes of <i>Shigella</i> , as well as enteroinvasive strains of <i>E. coli</i> , are disclosed.		

FOR THE PURPOSES OF INFORMATION ONLY

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AT	Austria	FI	Finland	MN	Mongolia
AU	Australia	FR	France	MR	Mauritania
BB	Barbados	GA	Gabon	MW	Malawi
BE	Belgium	GB	United Kingdom	NL	Netherlands
BF	Burkina Faso	GN	Guinea	NO	Norway
BG	Bulgaria	GR	Greece	NZ	New Zealand
BJ	Benin	HU	Hungary	PL	Poland
BR	Brazil	IE	Ireland	PT	Portugal
CA	Canada	IT	Italy	RO	Romania
CF	Central African Republic	JP	Japan	RU	Russian Federation
CG	Congo	KP	Democratic People's Republic of Korea	SD	Sudan
CH	Switzerland	KR	Republic of Korea	SE	Sweden
CI	Côte d'Ivoire	LI	Liechtenstein	SK	Slovak Republic
CM	Cameroon	LK	Sri Lanka	SN	Senegal
CS	Czechoslovakia	LU	Luxembourg	SU	Soviet Union
CZ	Czech Republic	MC	Monaco	TD	Chad
DE	Germany	MG	Madagascar	TG	Togo
DK	Denmark	MI	Mali	UA	Ukraine
ES	Spain			US	United States of America

NUCLEIC ACID PROBES FOR THE DETECTION OF SHIGELLA

Description

Background of the Invention

The genus *Shigella* includes four species (major
05 serogroups): *S. dysenteriae* (Grp. A), *S. flexneri*
(Grp. B), *S. boydii* (Grp. C) and *S. sonnei* (Grp. D)
as classified in Bergey's Manual for Systematic
Bacteriology (N.R. Krieg, ed., pp. 423-427 (1984)).
These serogroups are further subdivided into
10 serotypes (Table 1). The genera *Shigella* and
Escherichia are phylogenetically closely related.
Brenner and others have suggested that the two are
more correctly considered sibling species based on
DNA/DNA reassociation studies (D.J. Brenner, et al.,
15 International J. Systematic Bacteriology, 23:1-7
(1973)). These studies showed that *Shigella* species
are on average 80-89% related to *E. coli* at the DNA
level. Also, the degree of relatedness between
Shigella species is on average 80-89%. *Shigella*
20 *boydii* serotype 13 is atypical in that it is only
65% related to other *Shigella* serotypes and
Escherichia.

The genus *Shigella* is pathogenic in humans; it,
causes dysentery at levels of infection of 10 to 100
25 organisms. By contrast, the majority of *E. coli*

-2-

(9000 O:H serotypes) are not associated with diarrheal disease. Pathogenic E. coli serotypes are collectively referred to as Enterovirulent E. coli (EVEC) (J.R. Lupski, et al., J. Infectious Diseases, 157:1120-1123 (1988); M.M. Levine, J. Infectious Diseases, 155:377-389 (1987); M.A. Karmali, Clinical Microbiology Reviews, 2:15-38 (1989)). This group includes at least 5 subclasses of E. coli, each having a characteristic pathogenesis pathway resulting in diarrheal disease. The subclasses include Enterotoxigenic E. coli (ETEC), Verotoxin-Producing E. coli (VTEC), Enteropathogenic E. coli (EPEC), Enteroadherent E. coli (EAEC) and Enteroinvasive E. coli (EIEC). The VTEC include Enterohemorrhagic E. coli (EHEC) since these produce verotoxins.

The pathogenesis of Enteroinvasive E. coli is very similar to that of Shigella. In both, dysentery results from invasion of the colonic epithelial cells followed by intracellular multiplication which leads to bloody, mucous discharge with scanty diarrhea.

Thus, detection of Shigella and EIEC is important in various medical contexts. For example, the presence of either Shigella or EIEC in stool samples is indicative of gastroenteritis, and the ability to screen for their presence is useful in treating and controlling that disease. Detection of Shigella or EIEC in any possible transmission vehicle such as food is also important to avoid spread of gastroenteritis.

-3-

Currently, presence of Shigella in stool samples is detected by cultivating an appropriately prepared sample on microbiological media under conditions favorable for growth of those bacteria.

05 The resulting colonies are then examined for microbiological and biochemical characteristics, a process that typically takes at least three days and does not permit processing large numbers of samples. However, hospitals do not test for the presence of
10 EIEC in stool because of the difficulty of serotyping which is necessary to identify the EIEC among the numerous, non-pathogenic E. coli normally present in stool.

Summary of the Invention

15 The present invention relates to methods of detection of bacteria of the genus Shigella and/or Enteroinvasive E. coli (EIEC) by use of a set of nucleic acid probes. The invention further relates to a set of Shigella specific chromosomal sequences
20 and fragments, which were isolated by subtractive hybridization against non-Shigella DNA, and to probes derived from the Shigella specific fragments. Additionally, probes were derived from a sequence from the Shigella ompA gene. In particular, a
25 series of probes, each approximately 40 nucleotides in length, were designed having specificity for Shigella or for Shigella and Enteroinvasive E. coli, and having utility in non-isotopic test formats which require amplification to achieve high
30 sensitivity. In addition, specific hybridization

-4-

probe sets were developed which are capable of detecting substantially all clinically significant serotypes of Shigella, including S. sonnei, S. flexneri, S. boydii, and S. dysenteriae, as well as
05 enteroinvasive strains of E. coli.

Probes or probe sets of the present invention can be used in a number of hybridization formats for the detection of Shigella species and/or Enteroinvasive E. coli. For example, the presence
10 of one or more Shigella species and/or one or more species of EIEC in a sample can be determined by lysing the cells in the sample, contacting the sample with a DNA probe set under conditions
15 Shigella and/or EIEC DNA, capturing the hybrids formed between the probes and the sample DNA, and detecting the hybrid complexes by a suitable method as an indication of the presence of Shigella or EIEC in the sample.

20 Brief Description of the Figures

Figure 1 is a flow chart illustrating a strategy for the isolation of Shigella specific DNA sequences from "target" DNA and of Shigella specific fragments from a library of target DNA clones.

25 Figure 2 illustrates the nucleotide sequence (SEQ ID NO:1) of Shigella specific fragment NT6 and some flanking sequence, and the locations and sequences of probes 1500 (SEQ ID NO:14), 1501 (SEQ ID NO:15) and 1911 (SEQ ID NO:16), which are derived
30 from these sequences.

-5-

Figure 3 illustrates the nucleotide sequence of Shigella specific fragment NT11-2 (SEQ ID NO:2), and the locations and sequences of probes 1682 (SEQ ID NO:17), 1683 (SEQ ID NO:18), 1708 (SEQ ID NO: 19),
05 and 1709 (SEQ ID NO:20), derived from the fragment.

Figure 4 illustrates the nucleotide sequence of Shigella specific fragments NT14 (SEQ ID NO:4) and NT15 (SEQ ID NO:3), comprising a version of the Class III repeat isolated from S. sonnei, and the
10 location and sequence of probes 1864 (SEQ ID NO:22) and 437 (SEQ ID NO:21), derived from these fragments. The sequences of three versions of the Class III repeat isolated from E. coli (E.c. 1; SEQ ID NO:5; and E.c. 2; SEQ ID NO:6 and SEQ ID NO:7)
15 and S. flexneri (S.f.; SEQ ID NO:8 and SEQ ID NO: 9) are also shown.

Figure 5 illustrates the sequence of Shigella specific fragment NT18-1a (SEQ ID NO:10), and the location and sequence of probes 1712 (SEQ ID NO: 23)
20 and 1713 (SEQ ID NO: 24), derived from the fragment.

Figure 6 illustrates the sequence of Shigella specific fragment NT19-2 (SEQ ID NO: 11), and the location and sequence of probes 1684 (SEQ ID NO: 25) and 1685 (SEQ ID NO: 26), derived from the
25 fragment.

Figure 7 illustrates a portion of the sequence of the S. dysenteriae (S.d.) ompA gene (SEQ ID NO: 12), and the location and sequence of probes 1706 (SEQ ID NO: 27) and 1707 (SEQ ID NO: 28), derived
30 from the ompA gene sequence. The E. coli (E.c.)

-6-

ompA gene sequence (SEQ ID NO: 13) corresponding to the same region is shown for comparison.

Detailed Description of the Invention

Previous investigators interested in developing
05 specific probes for Shigella have targeted the
virulence plasmid. Both Shigella and EIEC harbor a
single copy virulence plasmid approximately 140 MD
(215 kilobasepairs) in size which is necessary for
invasion (T.L. Hale, Infection and Immunity,
10 40:340-350 (1983)). For example, U.S. Patent No.
4,816,389 (Sansonetti et al.) discloses a 27
kilobasepair (kb) region of the virulence plasmid
proven necessary for invasion. These investigators
have shown that the virulence plasmid, which exists
15 as one copy per bacterium, is unstable; regions of
the plasmid may be deleted. Shigella and EIEC
strains, which have been stored or passaged in the
laboratory, frequently are found to contain
virulence plasmids of reduced size (C. Sasakawa, et
20 al., Infection and Immunity, 51:470-475 (1986); A.T.
Maurelli, et al., Infection and Immunity,
43:397-401)). The 27 kb region from which the
probes of Sansonetti et al. have been derived has
been shown to be one of the unstable segments (P.K.
25 Wood, J. Clinical Microbiology, 24:498-500 (1986)).
Strains of Shigella or EIEC which do not contain a
portion of the 27 kb target region are not detected
by the probe and are incorrectly identified as
non-Shigella or non-EIEC.

-7-

Moreover, the 27 kb probe region contains insertion element 1 (IS1) which is ubiquitous among the Enterobacteriaceae (M. Venkatesan, et al., J. Clinical Microbiology, 26:261-266 (1988)). It also
05 contains at least one copy of insertion element 600 (IS600; S. Matsutani, et al., J. Molecular Biology, 196:445-455 (1987)), which occurs frequently in both pathogenic and non-pathogenic representatives of E. coli, as well as in Shigella and EIEC (unpublished
10 result). The presence of these broadly distributed insertion elements and the large size of the probes designed from the 27 kb region decrease the utility of these probes in non-isotopic test formats which require amplification in order to achieve high
15 sensitivity.

In contrast, the present invention relates to probes and probe sets which are (1) developed from Shigella specific fragments derived from chromosomal sequences of Shigella and (2) moderate in size, each
20 probe being approximately 40 bases in length. The increased stability of the chromosomal sequences detected by the probes compared to sequences of the virulence plasmid can result in increased reliability of detection. Furthermore, moderately
25 sized probes have utility in non-isotopic test formats which require amplification to achieve high sensitivity. Both the Shigella specific fragments and the probes derived from them are also useful as hybridization probes in other formats. Some of the
30 fragments derived from the Shigella chromosome are

-8-

also present in an episomal location on the invasion plasmid or another plasmid.

In one embodiment, the invention features a nucleic acid probe set consisting essentially of
05 nucleic acids with sequences that are:

- a. derived from the chromosomal sequence of representative bacteria of the species Shigella sonnei (ATCC 29930, designated type strain) and
10 Shigella flexneri type 2a (ATCC 29903, designated type strain) but are less than the entire chromosomal sequence of these bacteria;
- b. capable of hybridizing to DNA of members of the four known Shigella species and to Enteroinvasive E. coli (EIEC);
- 15 c. not capable of hybridizing or only weakly hybridizing to DNA of bacteria that are in neither the genus Shigella nor the group EIEC.

As used herein, a sequence fragment or oligonucleotide that is "derived from a chromosomal
20 sequence" is a natural, engineered or synthetic molecule having a sequence which is identical or complementary to a chromosomal sequence or is identical or complementary to a variant of the chromosomal sequence. A sequence fragment or
25 oligonucleotide which is identical or complementary to a variant of a selected chromosomal sequence (i.e., the variant differs in sequence from the

-9-

chromosomal sequence) is a homologue of the sequence or oligonucleotide which is identical or complementary, respectively, to the selected chromosomal sequence. This type of homologue will
05 have a nucleotide sequence substantially similar to a chromosomal sequence and will retain the desired function (will be able to hybridize to substantially the same nucleic acids as the sequence fragment or oligonucleotide which is identical or complementary
10 to that chromosomal sequence under similar hybridization conditions) of the fragment or oligonucleotide which is identical to the selected chromosomal sequence. A homologue may differ from the chromosomal sequence in sequence and/or may
15 contain modified nucleotides or nucleotide analogs (e.g., phosphorothioates, methylphosphonates). A homologue must be able to hybridize to the same nucleic acid as the sequence fragment or oligonucleotide which is identical or complementary
20 to that chromosomal sequence under similar hybridization conditions. It is well known to those skilled in the art that either strand of a double-stranded DNA sequence can serve as the target for a complementary probe. The complement of a
25 given probe is expected to have a substantially similar hybridization pattern, under similar hybridization conditions. The probes may be DNA or RNA or modified DNA or RNA.

Nucleic acid fragments or oligonucleotides
30 containing sequences derived from a chromosomal sequence, their homologues, and complements of all

-10-

of the foregoing may be used as probes. For example, the Shigella specific fragments, portions thereof and oligonucleotides from the Shigella specific fragments may be used as hybridization
05 probes.

In one embodiment, either of two probe sets of synthetically produced nucleic acid probes (each approximately 40 nucleotides long) will detect substantially all clinically significant serotypes
10 of Shigella, including S. sonnei, S. flexneri, S. boydii and S. dysenteriae. Clinically significant serotypes or isolates are those which are human pathogens. In addition, these probe sets recognize some or all enteroinvasive strains of E. coli,
15 exclusive of other enteric bacteria tested, with the exception of Escherichia fergusonii.

The invention further relates to methods of detecting Shigella species or EIEC in a sample. For example, one or more Shigella serotypes and/or EIEC
20 present in a sample can be detected by lysing the cells in the sample; contacting the sample with a nucleic acid probe or probes of the present invention under conditions that allow the probes to hybridize to Shigella and EIEC DNA in the sample,
25 thus forming hybrid nucleic acid complexes; isolating hybrid nucleic acid complexes formed between the probes and DNA in the sample, and detecting the hybrid nucleic acid complexes as an indication of Shigella or EIEC in the sample. For
30 example, clinical (e.g., stool), environmental (e.g. water), or food specimens may be subjected to such a

-11-

procedure to ascertain the presence of Shigella and/or EIEC species in the sample.

In preferred embodiments of the method, one or more pairs of probes are selected as a capture and
05 detector probe pair for use in a dual probe liquid hybridization format. These probes can be produced synthetically by chemical or enzymatic synthesis methods. They may be produced as part of a larger molecule. For example, the capture probe can be
10 tailed with 150-200 deoxyadenosine triphosphate (dATP) residues using the enzyme terminal deoxynucleotidyl transferase. The detector probe can be incorporated into an amplification/detection system, such as a biotin-streptavidin-alkaline
15 phosphatase system. Both the capture and detector probes are then allowed to hybridize to the target nucleic acid in a background of competitor nucleic acid from the sample. The hybrid products can be captured out of the mixture by magnetic beads
20 complexed with tails of deoxythymidine monophosphate generally 14 residues longer. The captured hybrids (i.e., those affixed to or captured on the magnetic beads) are detected by means of the amplification/detection system selected.

25 Isolation of Shigella-Specific Fragments

Shigella-specific genomic sequences were isolated by subtractive hybridization using biotin-streptavidin agarose affinity column chromatography methods described by Langer et al.
30 (Langer, P.R. et al., Proc. Natl. Acad. Sci. USA,

-12-

78:6633-6637 (1981) and Welcher et al. (Welcher, A.A. et al., Nucleic Acids Res., 14: 10027-10044 (1986)), each of which is incorporated herein by reference. The method is outlined in Figure 1.

05 For example, a mixture of competing DNAs, such as the complex DNA competitor mix described in Figure 1 can be used. In general, a complex DNA competitor mix will contain an excess of DNA (relative to the target DNA) from one or more
10 Enterohemorrhagic E. coli isolates, one or more Enterotoxigenic E. coli isolates, and one or more non-pathogenic E. coli isolates. The mix may further contain specific sequences already recovered or sequences designed to eliminate non-specific
15 sequences. For example a plasmid containing a 17 kb region of the Shigella virulence plasmid such as pHS4033, Class III repeat DNA, or M13mp18 RF DNA may be included.

The particular complex DNA competitor mix used
20 contained a 6-fold excess by weight of DNA (relative to the target DNA) from each of Enterohemorrhagic E. coli isolate IG 3040, Enterotoxigenic E. coli isolate IG 3060, and non-pathogenic E. coli (YMC). The mix further contained DNA, in amounts equal to
25 the target DNA by weight, from each of the following: a pBR322 clone containing a 17 kb region of the Shigella virulence plasmid (plasmid pHS4033, Boileau, C.R. et al., J. Clin. Microbiol. 20(5): 959-961 (1984)), M13mp18 RF DNA, and the sequence of
30 the Class III repeat (Class IIIR-IG900) cloned into pBR322. (The 1.3 kb Class III repeat and adjacent

-13-

chromosomal DNA was cloned into pBR322. The length of the insert was 3.5 kb.) The complete mixture was labelled with biotin-11-dUTP by the nick-translation method. The "target DNA" from which specific
05 sequences were identified was isolated from a single *Shigella* species, digested with restriction endonuclease *Sau3A*, and end-labelled with ^{32}P in an end-filling reaction with DNA polymerase I.

The two DNA pools were combined such that the
10 competing DNAs were at a 20-40 fold molar excess relative to the *Shigella* DNA. The mixture was denatured and hybridized in liquid at low stringency overnight. The hybridization buffer contained 0.75 M NaCl, 50 mM NaPO_4 , 1 mM EDTA, 0.05% SDS and 40%
15 formamide. The hybridization mixture then was passed over the streptavidin agarose column. *Shigella* DNA sequences that were sufficiently complementary to the competing DNAs to form hybrids under the conditions used were retained on the
20 streptavidin agarose affinity column by virtue of the biotin incorporated into the competitor DNA. In contrast, sequences that were unable to form hybrids under those conditions were enriched in the nucleic acid fraction passing through the column. The
25 latter sequences contain *Shigella*-specific sequences.

A small aliquot of the DNA enriched for *Shigella*-specific sequences (previously ^{32}P -labelled) was used to probe *Shigella* and *E. coli*
30 DNAs of interest (0.1 μg of each DNA was spotted in a 3 μl volume on nitrocellulose). Hybridization

-14-

conditions were those described below in Example II for nick-translated fragments. Cross-hybridization to the competitor E. coli provided an indication that further enrichment of the Shigella "target" DNA was required. Typically, four cycles of competition hybridization versus the biotinylated E. coli competitor DNA were necessary to eliminate the cross-hybridization. This was accomplished by repeating the competition hybridization/affinity capture cycle, using the Shigella DNA which passed through the column during the previous cycle, as the starting material for the next cycle. In this way, the labeled Shigella target DNA was progressively enriched for Shigella-specific sequences in each cycle.

The nucleic acid which was enriched for Shigella-specific sequences was then used to probe a Sau3A library of the same Shigella isolate used as the "target DNA" in the subtractive hybridization; the library was constructed in the plasmid vector, pUC18. Inserts (fragments) from positive clones were purified from the vector, labelled with ^{32}P by nick-translation and used to probe mini-cyto-dot panels of inclusivity and exclusivity organisms as described in Examples I and II. A probe shows inclusivity toward an organism if DNA from the organism hybridizes to the probe, and shows exclusivity toward an organism if the probe does not hybridize (or if hybridization is barely detectable) to DNA from that organism under the particular hybridization conditions used.

-15-

In most cases, hybridization to inclusivity organisms was much stronger than hybridization to competitors or exclusivity organisms. When necessary, subcloning was done to remove the
05 sequences which cross-hybridized to competitors. At this point, the fragments were labelled with ^{32}P by nick-translation and used to probe full inclusivity and exclusivity cyto-dot panels.

The inclusivity panels consisted of bacteria
10 representing all known Shigella serotypes as well as Enteroinvasive E. coli which exhibit the same pathogenesis as Shigella. These organisms are listed in Tables 2 through 6. The exclusivity panels consisted of non-pathogenic E. coli,
15 Enterotoxigenic E. coli, Enteropathogenic E. coli, Enterohemorrhagic E. coli, other Escherichia species and gram negative Enterobacteriaceae. The exclusivity organisms are listed in Tables 6 and 7a/7b.

20 In the Tables, under the conditions used, (-) indicates no signal, (+/-) or (-/+) was barely detectable, (+) indicates a weak but reproducible and readily detectable signal, (++) indicates a moderate signal, (+++) indicates a strong signal,
25 and (+++++) indicates a signal comparable to the positive control (DNA from the organism from which the fragment was isolated, or a known sequence identical to the probe being tested). The genomic (chromosomal) DNA fragments identified by the
30 subtractive hybridization and refinement protocols

-16-

described above are referred to as Shigella specific fragments.

1. Shigella specific fragments NT19-2 (SEQ ID NO: 11) and NT18-1a (SEQ ID NO: 10) were isolated from a library of Shigella flexneri (ATCC 29903) genomic clones by probing the library with Shigella specific sequences isolated from S. flexneri (ATCC 22903) ³²P-labeled "target" DNA, using the "complex DNA competitor mix" described above in the subtractive hybridization steps.
2. Shigella specific fragments NT 6 (see SEQ ID NO:1), NT11-2 (SEQ ID NO:2), NT14 (SEQ ID NO: 4) and NT15 (SEQ ID NO:3) were isolated from a library of Shigella sonnei (ATCC 29930) genomic clones by probing the library with Shigella specific sequences isolated from S. sonnei (ATCC 29930) ³²P-labeled "target" DNA and using a single non-pathogenic E. coli (YMC) and 1 µg of pBR322 vector DNA as source of competitor DNA in the subtractive hybridization. In this case, avidin agarose rather than streptavidin agarose was used in the affinity column.
- 25 The hybridization results against inclusivity and exclusivity organisms for these Shigella-specific fragments are recorded in Tables 2-7. In addition, a set of ompA probes were developed from published

-17-

data (G. Braun, et al., Nucleic Acids Research,
10:2367-2378 (1982)) by comparison of the outer
membrane protein gene sequence (ompA) of Shigella
dysenteriae (SEQ ID NO:12) with that of E. coli (SEQ
05 ID NO:13). The regions of ompA gene sequence which
appeared most different between the two sequences
were selected for the development of test probes.
These were synthesized and assayed by hybridization
analysis as described above.

10 Isolation of Shigella-Specific Oligonucleotides

The Shigella-specific fragments described above
displayed the most marked specificity for
inclusivity organisms of the Shigella DNA fragments
analyzed. These fragments were sequenced and
15 oligonucleotide probes useful as capture and
detector probes were designed from these sequences.
Oligonucleotide capture and detector probes were
also designed from a fragment of the ompA sequence
by comparative sequence analysis. Following
20 synthesis, the oligonucleotides were end-labelled
with ^{32}P and tested by cyto-dot hybridization
analysis as described in Example I to ensure that
they exhibited the desired inclusivity and
exclusivity behavior or pattern. In several cases,
25 an additional exclusivity panel was tested at this
point (Table 8). This panel consisted of 4 μg DNA
dots of gram positive and gram negative
bacteria--including aerobic and anaerobic
representatives commonly found in stool. The DNA
30 was isolated by a protocol which makes use of glass

-18-

beads to physically disrupt the cell wall of bacteria. Each DNA was spotted on nitro-cellulose filters in a 3 μ l volume; the DNAs were denatured, neutralized and fixed as described in Example I for
05 the preparation of cyto-dot panels.

Description of Probes and Hybridization Behavior
With Respect to Inclusivity and Exclusivity

Oligonucleotide probes were derived from the Shigella-specific fragments identified by affinity
10 chromatography and from the ompA sequence. The sequence of each oligonucleotide probe and the Shigella-specific DNA fragment that each was derived from are listed in Table 9. The inclusivity and exclusivity hybridization behavior of the clones and
15 subclones, and the oligonucleotides designed from these sequences are described below.

Fragment NT6 and Probes 1500, 1501, and 1911

NT6 (see SEQ ID NO:1) is a 124 bp Sau3A Shigella specific fragment. The sequence is
20 repeated 6 times in Shigella sonnei (ATCC 29930) chromosomal DNA. It is also found in one or two copies on the virulence plasmid of other Shigella isolates. The entire fragment was sequenced (Figure 2). The first 124 bp of Figure 2 are from NT6 (see
25 SEQ ID NO:1). The Sau3A site at the 3'-end of the NT-6 sequence is indicated. In this and other Figures, IUPAC conventions for referring to nucleotides and sequence ambiguities are used:

-19-

Unambiguous bases: A, C, G, and T (or U)

Two possible bases: M (A or C)

R (A or G)

W (A or T)

05

S (C or G)

Y (C or T)

K (G or T)

Three possible bases: B (C, G or T)

D (A, G or T)

10

H (A, C or T)

V (A, C or G)

Four possible bases: N (A, C, G or T)

Where there is a weak band on the sequencing gel, but no other band, the base is indicated by a lower case letter. Regions of ambiguous nucleotide order due to band compression are enclosed in parentheses.

Two oligonucleotides derived from NT6, each 35 bases long (1500, SEQ ID NO:14) and 1501, SEQ ID NO: 15), were designed and can be used as capture and detection probes (Table 9). A third probe, 40 bases long (1911, SEQ ID NO:16), was designed from NT6 and from additional sequence adjacent to the Sau3A NT-6 fragment (38 bp; see SEQ ID NO:1). This additional sequence was obtained from a clone isolated from an S. sonnei library using the NT-6 124 bp fragment as a probe. A sequencing primer internal to NT-6 was used for sequencing. Oligonucleotide 1911 was

-20-

designed and tested. The hybridization pattern of 1911 was identical to that of NT6, suggesting it is part of the repeated element.

The individual oligonucleotides were tested for
05 hybridization to inclusivity and exclusivity
organisms. In cases where all three probes
displayed the same hybridization pattern, the
results are not listed separately. The individual
oligonucleotides, like the parent fragment,
10 hybridized (using a + signal as the lower limit of
hybridization) to all S. sonnei isolates tested,
many S. dysenteriae and S. boydii, S. flexneri
serotypes, including S. flexneri type 6 (Tables 2
through 5 and summary Table 10). The
15 oligonucleotide probes derived from NT6 hybridized
to all Enteroinvasive E. coli, but not to other
classes of pathogenic E. coli, non-pathogenic E.
coli or other organisms commonly found in stool
(Tables 6-8).

20 Fragment NT11-2 and Probes 1682, 1683, 1708 and 1709

NT11-2 (SEQ ID NO:2) is a 796 bp Hha I subclone
of an original Sau3A fragment (NT11) which was 3.5
kb in length. Fragment NT11-2 has been sequenced
(Figure 3; SEQ ID NO:2). Two sets of
25 oligonucleotide probes, useful as a capture/detector
probe pair were designed from opposite ends of the
fragment. Probes 1682 (SEQ ID NO:17) and 1683 (SEQ
ID NO:18) are 41 nucleotides long (Table 9). Probes
1708 (SEQ ID NO:19) and 1709 (SEQ ID NO:20) are 35
30 and 36 nucleotides long, respectively. When used as

-21-

capture/detector probes, each capture probe can be synthesized with an additional three nucleotides (TAT) at the 3' end to ensure efficient tailing by the enzyme terminal deoxynucleotidyl transferase.

05 Tables 2 through 5 and summary Table 10 show that the two sets of oligonucleotides have different hybridization patterns. The cumulative hybridization pattern using the four probes together is the same as that of the parental fragment, NT

10 11-2, which hybridizes to all S. sonnei isolates, and to some S. dysenteriae, S. flexneri and S. boydii isolates. However, when used as capture/detector probe pairs (i.e., 1682 paired with 1683 or 1708 paired with 1709), certain serotypes

15 would not be detected because one partner of each probe set does not hybridize to the isolate (e.g., S. dysenteriae types 9 and 10, S. boydii type 17). One of the four Enteroinvasive E. coli is detected strongly with 1708/1709 and weakly with 1682/1683.

20 Two pathogenic E. coli isolates would be detected by 1708/1709 while only one of these isolates would be detected with 1682/1683. Apart from these pathogenic isolates, the oligonucleotides do not cross-hybridize to other non-pathogenic E. coli or

25 other organisms commonly found in stool (Tables 6 and 7a).

Fragments NT14 and NT15, and Probes 437 and 1864

NT14 (SEQ ID NO:4) and NT15 (SEQ ID NO:3) are Sau3A fragments which are approximately 800 bp and

30 600 bp in size, respectively. The two fragments

-22-

have been sequenced and are portions of a highly repeated element which is 1.3 kb in length. One representative repeat each from S. sonnei (ATCC 29930) and S. flexneri (ATCC 29903) and two repeats
05 from E. coli (IG900) were cloned and sequenced (Figure 4; S. sonnei repeat, SEQ ID NO: 3 and SEQ ID NO:4; E. coli repeat 1, SEQ ID NO:5; E. coli repeat 2, SEQ ID NO:6 and SEQ ID NO:7; S. flexneri repeat, SEQ ID NO:8 and SEQ ID NO:9). The repeat is highly
10 conserved and has characteristics of a transposable element. Over 20 copies of the repeat sequence are present in the chromosome and virulence plasmid of Shigella. The repeat occurs in 1 to 3 copies in some E. coli competitors, but not in other bacterial
15 species.

There are only a few differences between the E. coli and S. sonnei sequences shown in Figure 4. A 17 base oligonucleotide probe (probe 1864, SEQ ID NO: 22) was designed such that a single mismatch is
20 located eight bases from either end of the probe (Table 9). This probe hybridizes strongly to the majority of Shigella and EIEC tested and does not cross-hybridize to competitor bacteria (Tables 2-8 and summary Table 10). A companion detector probe
25 can be designed within the boundaries of the Class III repeat on either side of the specific capture probe, 1864. One such example is the complement of probe 437 (SEQ ID NO:21) (49 bases) for which inclusivity and exclusivity is listed in Tables 2
30 through 8. This probe hybridizes strongly to all serotypes of Shigella except S. boydii serotype 13.

-23-

Probe 437 hybridizes weakly to a number of non-pathogenic and pathogenic strains of E. coli, but does not hybridize to other Enterobacteriaceae tested. The hybridization signal with certain E. coli is due to the low copy number of the Class III repeat in this genus versus the high copy number of the repeat in the genus Shigella.

Two extra bases incorporated into the sequence of probe 437 as a result of sequencing errors in the S. sonnei sequence (a G and a T) have proved useful in decreasing the signal for E. coli isolates relative to Shigella strains. The 437 probe does not hybridize as well as expected to the positive control for S. sonnei or S. flexneri, suggesting that the two nucleotides in question are not present in the S. sonnei sequence. This observation also is likely to be related to copy number differences between the two genera.

Probe 1864 (SEQ ID NO:22) hybridized to all isolates of S. dysenteriae tested except for one isolate of serotype 1 (IG 826). However, serology on this isolate was not confirmed.

Fragment NT18-1a and Probes 1712 and 1713

NT 18-1a (SEQ ID NO:10) was subcloned from the original Sau3A fragment, NT18, in two steps. A PstI/SacI double digest of NT18 (1500 bp) yielded fragment NT18-1 (1250 bp), which was then restricted with HaeIII to generate NT18-1a (630 bp). Sequences related to fragment NT18 are also known to occur on a small multicopy plasmid which is distinct from the

-24-

215 kbp invasion plasmid. The sequence of NT18-1a is shown in Figure 5 (SEQ ID NO:10). Oligonucleotide probes 1712 (SEQ ID NO:23) and 1713 (SEQ ID NO:24), suitable as capture/detection probes, both 37 bases long, were designed from the sequence (Figure 5 and Table 9). The hybridization pattern (isolates detected) of the oligonucleotide probes was identical to the parental *Shigella* specific fragment NT18-1a with the exception of one *S. flexneri* isolate (Tables 2-5, 7b). This strain (IG 711) was detected by oligonucleotide probe 1713, but not by probe 1712. In a liquid hybridization assay when the two probes would be used as a capture/detection probe pair this organism would not be detected. Under the conditions used, the probes hybridize to 6/8 type 1 *S. dysenteriae*, to all *S. flexneri* isolates with the exceptions of three type 6 isolates, the IG711 isolate mentioned above, IG872, IG741, and IG709 (Tables 2 through 5 and summary Table 10). The probes do not detect Enteroinvasive *E. coli*, but cross-hybridize to one pathogenic *E. coli* under the conditions used. They do not cross-hybridize to non-pathogenic *E. coli* or other bacteria commonly found in stool (Tables 6 through 8).

Fragment NT19-2 and Probes 1684 and 1685

Fragment NT19-2 (388 bp; SEQ ID NO:11) is an *RsaI* subclone of the original *Sau3A* fragment which was 1070 bp in length. NT19-2 was sequenced (Figure 6; SEQ ID NO:11) and oligonucleotide probes 1684

-25-

(SEQ ID NO:25) and 1685 (SEQ ID NO:26), each 35 nucleotides long, were designed (Table 9). These probes are suitable as capture/detection probes. The hybridization patterns or spectrum of isolates
05 detected by the individual oligonucleotides and the parental fragment are identical. Hybridization to some S. boydii, some S. sonnei and all S. flexneri except type 6 was observed (Tables 2-5, summary Table 10). The probes hybridize to one out of five
10 Enteroinvasive E. coli, and do not cross-hybridize to other pathogenic E. coli, non-pathogenic E. coli or other bacteria commonly found in stool (Tables 6-8).

ompA Fragment and probes 1706 and 1707

15 Oligonucleotides 1706 (SEQ ID NO:27) and 1707 (SEQ ID NO:28) were designed from the published sequence of the outer membrane protein gene (ompA) of Shigella dysenteriae. Figure 7 shows the S. dysenteriae ompA gene sequence from nucleotide
20 position 893 through 1076, according to the numbering of Braun et al. (Nucl. Acids Res. 10(7): 2367-2378 (1982); SEQ ID NO:12). This region contains significant differences between the E. coli and S. dysenteriae ompA coding sequences. The
25 sequence of the corresponding region of the E. coli ompA gene is shown for comparison (SEQ ID NO:13), and the positions of probes 1706 (SEQ ID NO:27) and 1707 (SEQ ID NO:28) are indicated in Figure 7.

Both oligonucleotides are 35 bases long (Table
30 9). Probe 1706 has 7 differences between the E.

-26-

coli and the S. dysenteriae sequence. The region from which probe 1707 was designed is 15 bases shorter in E. coli. Additionally, there are numerous differences in that probe site between E. coli and S. dysenteriae. These probes hybridize with S. dysenteriae types 1 and 2 and many S. boydii serotypes (Tables 2-5, summary Table 10). When probe 1707 is used as the specific capture probe and 1706 is used as the detector probe in liquid hybridization, no hybridization is anticipated to Enteroinvasive E. coli, other pathogenic E. coli, non-pathogenic E. coli and other bacteria commonly found in stool with the exception of Escherichia fergusonii (Tables 6-8).

15 Description of Probe Sets

Probe Set I

A desired inclusivity/exclusivity pattern may be achieved by use of various combinations of probes. One possible strategy involves pooling one or more combinations of probe pairs to make a probe set. For example, capture/detection oligonucleotides designed from fragments NT 6 (SEQ ID NO:1), NT 19-2 (SEQ ID NO:11) and the ompA gene (SEQ ID NO:12) may be used together as a probe set or combination for detection of substantially all clinically significant Shigella serotypes. One possible probe set comprises three capture/detection probe pairs, including probe pairs 1684/1685, 1707/1706, and a pair selected from probes

-27-

1500/1911/1501. This substantially inclusive probe set detects all *Shigella* serotypes except for *S. dysenteriae* type 10 and *S. boydii* type 13 (using a + signal as a lower limit for hybridization). In addition, this probe set does not cross-hybridize under the conditions used to any of the competitors tested except for *Escherichia fergusonii*.

The two *Shigella* serotypes that are not detected by probe set I under these conditions, are rarely isolated in the United States, as indicated by records of the Center of Disease Control. For example, out of a total of 167,915 cases of *Shigella* infection reported by the Center for Disease Control for the years 1976 through 1987, only two cases were identified as *S. dysenteriae* type 10 and three cases as *S. boydii* type 13.

In a capture/detection assay format, the more specific oligonucleotide of a capture/detector pair is preferred as the capture probe. In the case of the 1684/1685 probe set, either oligonucleotide may serve as the capture probe with equivalent hybridization results. However, in the case of the 1707/1706 probe set, oligonucleotide 1707 is preferred as the capture probe since it does not cross-hybridize to competitors. (Probe 1706 has (+/-) hybridization signal with certain competitors and a strong signal with an *E. blattae* isolate (Tables 6, 7b)).

In the case of the 1500/1911/1501 oligonucleotide probes, it is best to use either 1500 or 1911 as the capture probe. Probe 1501 has a

-28-

(+/-) hybridization signal with certain competitors (Table 8), but may be used as a detector probe with no adverse effects.

Probe Set II

05 An alternative probe set combines
capture/detection oligonucleotide probe pairs from
the Class III repeat (fragments 14 and 15; SEQ ID
NO:4 and SEQ ID NO:3, respectively) and the ompA
gene (SEQ ID NO:12). This substantially inclusive
10 probe set 1707/1706 and 1864/437-complement)
hybridizes to all *Shigella* except *S. boydii* type 13
and cross-hybridizes weakly to *Escherichia*
fergusonii.

 In the case of the capture/detector pair
15 1864/437-complement, it is best to use
oligonucleotide 1864 as the capture probe since it
is the more specific probe of the pair.
Oligonucleotide probes designed from sequences to
the left or right (in the 5' or 3' direction and
20 from the same strand) of the sequences from which
probe 1864 was derived may also serve as detector
probes. The complement of probe 437 is one such
example, and is expected to substantially retain the
hybridization pattern of probe 437.

25 Use of probe combination II requires only four
oligonucleotides instead of six, yet gives the
desired inclusivity and exclusivity (substantially
inclusive, in this case). In addition, the target
of probe 1864 is present in multiple copies (20-30
30 copies), and therefore, allows for increased

-29-

sensitivity. However, carefully controlled hybridization conditions are necessary to maintain exclusivity with probe set II, since the specificity of probe 1864 depends on a single mismatch to
05 differentiate between *Shigella* and *E. coli* which harbor the repetitive element (Class III repeat).

Table 12 lists the number of isolates expected to be detected by probe set I and probe set II in a capture/detection format using a (++) hybridization
10 signal as the cut-off for detection. The results for probe sets I and II in a capture/detection format are expected to be the same when a (+) signal is used as the cut-off for detection.

Additional Probes

15 Other probes (double- or single-stranded nucleic acid fragments or oligonucleotides), probe sets or combinations may be derived from the *Shigella* specific fragments. These fragments or oligonucleotides (probes) "derived from the sequence
20 of *Shigella* specific fragments", comprise nucleic acid sequences which are identical or complementary to a portion of the sequence of the *Shigella* specific fragments (and therefore to the *Shigella* chromosome). In some cases, only a portion of the
25 probe may be identical to the sequence of the original *Shigella* specific fragment. Portions of a probe which are identical or complementary to the sequence of a *Shigella* specific fragment can be noncontiguous in the probe.

-30-

The preferred probes will retain features of the inclusivity, and if desired, the exclusivity of the Shigella specific fragments. A probe derived from a Shigella specific fragment which

- 05 "substantially retains the inclusivity behavior of" a selected Shigella specific fragment hybridizes, under the same conditions, to at least one isolate of 90% or more of the (typed) serotypes to which the original fragment hybridizes. An original Shigella
- 10 specific fragment which hybridizes to at least one isolate of one of the 35 serotypes listed in Table 1 is said to hybridize to that serotype. A probe which "moderately retains the inclusivity behavior of" a selected Shigella specific fragment
- 15 hybridizes, under the same conditions, to at least one isolate of 83% or more, but less than 90%, of the serotypes to which the original fragment hybridizes. A probe which "partially retains the inclusivity behavior of" a selected Shigella
- 20 specific fragment hybridizes, under the same conditions, to at least one isolate of 50% or more, but less than 83%, of the serotypes to which the original fragment hybridizes.

- Exclusivity was determined using two sets of
- 25 exclusivity organisms. The exclusivity organisms screened included 152 non-EIEC strains listed in Tables 6, 7A and 7B, and defined here as non-EIEC Enterobacteriaceae exclusivity organisms. In addition, the 91 strains listed in Table 8, comprise
- 30 a second set of exclusivity organisms defined here as exclusivity organisms commonly found in stool.

-31-

A probe derived from a *Shigella* specific fragment which has "improved" exclusivity behavior for a given set of exclusivity organisms (e.g., non-EIEC Enterobacteriaceae) is one for which all of

05 the exclusivity organisms of that set have been screened in the dot blot format, and which detects (hybridizes with a signal of (+) or better) fewer exclusivity organisms of that set than the *Shigella* specific fragment from which it is derived, under

10 the same hybridization conditions. A probe derived from a *Shigella* specific fragment which, for a given set of exclusivity organisms, "substantially retains" the exclusivity behavior of the fragment from which it is derived, is one for which 90% or

15 more of the exclusivity organisms of that set have been screened in the dot blot format, and which has substantially the same or identical exclusivity behavior under the same hybridization conditions. In particular, a probe which will detect no more

20 than 13 strains of a set of exclusivity organisms which are not detected by the fragment from which it is derived, and which may or may not detect the exclusivity organisms which are detected by the original fragment, is defined as one which

25 "substantially retains" the exclusivity behavior of that fragment. It will be appreciated that a probe for which exclusivity has been determined for 100% of a given set of exclusivity organisms, but which detects the same number of exclusivity strains or

30 more (but not more than 13 additional strains) of

-32-

the exclusivity organisms as the original fragment falls in this latter category.

Furthermore, a probe for which exclusivity has not been determined for 100% of the organisms may be shown to have "improved" exclusivity behavior. For example, probes 1684 (SEQ ID NO:25) and 1685 (SEQ ID NO:26), for which the exclusivity for 4 non-EIEC Enterobacteriaceae has not been determined, do not detect two strains in Table 6 which are detected by NT19-2 (SEQ ID NO:11). If it is determined that these probes do not detect three or four of the strains not tested, then they will have improved exclusivity behavior, although they are presently classified as substantially retaining the exclusivity of NT19-2. Thus, the two classifications are not mutually exclusive.

Homologues of fragments and oligonucleotides derived from Shigella specific fragments, which hybridize to substantially the same serotypes as the fragments and oligonucleotides derived from Shigella specific fragments under the same hybridization conditions, can also be used. Homologues of a sequence fragment or oligonucleotide derived from a Shigella specific sequence will be identical or complementary to all or part of a variant of the sequence of a Shigella specific fragment.

For example, oligonucleotide probes, typically from about 10 nucleotides in length up to about 50 nucleotides in length, comprising sequences identical to a portion of a Shigella specific fragment can be designed. However, an

-33-

oligonucleotide probe may be longer than 50 nucleotides. Larger fragments comprising a sequence identical to a portion of a full length fragment can be prepared by restriction digestion of an isolated
05 clone, exonuclease digestion, by the polymerase chain reaction using selected primers, or other suitable methods, for example.

The additional probes derived from the Shigella specific fragments, complements or homologues
10 thereof, can be screened using a dot blot format (such as the cyto-dot or DNA dot blot formats of the Examples). These additional fragments or oligonucleotide probes can be used alone or in various probe pairs or probe combinations, or in
15 addition to a selected probe, probe pair or combination from Table 9. The additional probes can also be used as alternatives to the probes listed in Table 9. For example, another probe derived from Shigella specific fragment NT11-2 (SEQ ID NO:2)
20 could be used together with probe 1682 (SEQ ID NO:17) in place of probe 1683 (SEQ ID NO:18). For use in a capture/detection format, the probe would be derived from the same strand of the Shigella specific fragments as probe 1682. A probe which
25 substantially retained the inclusivity pattern of NT11-2 could be selected, for example.

Furthermore, probes of the present invention can be used in combination with other probes for Shigella, enteroinvasive E. coli, or other organisms (e.g.,
30 Salmonella, Campylobacter, etc.).

-34-

It will be appreciated that the signal for the recommended probe sets may be increased by using additional probes. Additional probes may be selected from the probes listed in Table 9,
05 oligonucleotide probes derived from the large Shigella specific fragments disclosed, the homologues and complements of any of the foregoing, or other suitable probes. Although some probes are preferred as detection probes in a capture-detection
10 probe format due to hybridization with exclusivity organisms, each probe may be used as either a capture or detection probe.

Inclusivity and Exclusivity Patterns

Different inclusivity and exclusivity patterns
15 can be obtained using selected combinations of probes. Furthermore, inclusivity and exclusivity behavior may be modulated by hybridization conditions, and/or by taking a specific level of hybridization as the cut-off. For example, in Table
20 10, a (+) signal, which is a weak but reproducible and readily detectable signal, is used as the cut-off for inclusivity or detection in the dot blot format (cyto-dot or DNA dot format).

In Table 11, a (++) cut-off is used. In Table
25 11, the number of isolates of each serotype or untyped isolate to which the probes NT6, NT11-2, NT18-1a, NT19-2, hybridized with a signal of at least (++) is indicated. In addition, the expected number of isolates of each serotype or untyped
30 isolates to which probe pairs selected from

-35-

1911/1500/1501, 1682/1683, 1708/1709, 1712/1713,
1684/1685, or 1706/1707 are expected to hybridize
with a signal of at least (++) in a
capture/detection format are also indicated. (ND
05 indicates that hybridization was not determined for a
particular isolate.)

Unless indicated otherwise herein, a probe
which "detects" or for "the detection" of an isolate
or serotype is one which gives at least a (+) signal
10 under the hybridization conditions used with the
isolate or with an isolate of a specific serotype.
An individual probe (fragments or oligonucleotides),
probe pair or probe set (a combination of probes
and/or probe pairs) which, under the conditions used
15 in the dot blot format, detects (hybridizes with a
signal of at least +) at least one isolate of 90% or
more of the serotypes listed in Table 1 is defined
as a "substantially inclusive" probe, probe pair, or
probe set. Similarly, a probe, probe pair, or probe
20 set which, under the conditions used in the dot blot
format, detects at least one isolate of 83% or more,
but less than 90%, of the serotypes listed in Table
1 is a moderately inclusive probe, probe pair or
probe set. A probe, probe pair or probe set which,
25 under the conditions used in the dot blot format,
detects at least one isolate of 50% or more, but
less than 83%, of the serotypes listed in Table 1 is
a partially inclusive probe, probe pair or probe
set.

30 For example, probe sets I and II described
above are substantially inclusive probe sets. Based

-36-

on the dot blot data (see summary Table 10), these combinations are expected to detect at least one member of 33 and 34 out of 35 serotypes, respectively in a capture/detection format (94.2% and 97.1%). In fact, each of these probe sets is expected to detect every isolate tested of the serotypes detected.

Individual probes such as oligonucleotides 1900, 1500 or 1501, which detect at least one member of 60% of the serotypes listed in Table 1, or fragment NT-6, which detects approximately 68% of the serotypes listed in Table 1 (see Table 10), would be considered partially inclusive probes. Probe 1864 is a substantially inclusive probe. Partially, moderately and substantially inclusive probes may be combined with each other or with other probes into appropriate pairs or sets to achieve a desired inclusivity pattern.

As stated above, the exclusivity organisms screened include the 152 non-EIEC strains listed in Tables 6, 7A and 7B, and defined here as non-EIEC Enterobacteriaceae exclusivity organisms. In addition, the 91 strains listed in Table 8, define a second set of exclusivity organisms defined here as exclusivity organisms commonly found in stool. An individual probe (fragments or oligonucleotides), probe pair or probe set (a combination of probes and/or probe pairs) which, under the conditions used in the dot blot format, does not detect (hybridizes with a signal less than +) any of the 152 non-EIEC strains listed in Tables 6, 7A and 7B is defined as

-37-

an "exclusive" probe, probe pair or probe set with regard to the non-EIEC Enterobacteriaceae exclusivity organisms. A probe, probe pair, or probe set which, under the conditions used in the
05 Tables, detects (hybridizes with a signal of (+) or better) 10% or less of the 152 non-EIEC strains listed in Tables 6 and 7A and B is defined as a "substantially exclusive" probe, probe pair or probe set with regard to these E. coli and
10 Enterobacteriaceae exclusivity organisms, while a probe which detects 20% or less of the 152 non-EIEC strains listed in Tables 6 and 7A and B is defined as "moderately exclusive" of the non-EIEC Enterobacteriaceae exclusivity organisms. A probe
15 which is "exclusive" also meets the criteria for a moderately or substantially exclusive probe.

A probe, probe pair or probe set which, under the conditions used in the dot blot format, does not detect (hybridizes with a signal less than +) any of
20 the 91 strains commonly found in stool and listed in Table 8 is defined as an "exclusive" probe, probe pair or probe set with regard to the exclusivity organisms commonly found in stool. A probe, probe pair, or probe set which, under the conditions used
25 in the Tables, detects (hybridizes with a signal of (+) or better) to 10% or less of the 91 strains listed in Table 8 is defined as "substantially exclusive" of the exclusivity organisms commonly found in stool. Thus, a probe that is exclusive of
30 the organisms in Table 8 is also substantially exclusive of the same organisms.

-38-

For example, probes 1500, 1501 or 1911 are exclusive of the non-EIEC Enterobacteriaceae of Tables 6 and 7, as well as of the strains commonly found in stool. Probes 1684 and 1685 did not detect
05 any of the non-EIEC tested; however, 4 organisms (2.6%) were not tested. Thus, 1684 and 1685 are known to be substantially exclusive. If none of these 4 organisms is detected, then these probes would be exclusive. Probe set I, when used in a
10 capture/detection format, is expected to detect 2 of the strains from Table 7 (Escherichia fergusonii). However, four strains from Table 7 were not tested (ND). Therefore, probe set I detects at most 3.9% (6/152) of the non-EIEC Enterobacteriaceae, and is
15 thus "substantially exclusive" of the non-EIEC Enterobacteriaceae exclusivity organisms. Probe set I is also exclusive of the exclusivity strains commonly found in stool (listed in Table 8).

Probe 437 detects about 9.9% (15) of the
20 non-EIEC organisms; however, 10/152 of the non-EIEC Enterobacteriaceae were not screened. Thus, the exclusivity of this probe for the non-EIEC Enterobacteriaceae is between 16.4% (moderately exclusive) and 9.9% (substantially exclusive).
25 However, in a capture/detection format, probe set II, which may use the complement of 437, is substantially exclusive of the non-EIEC Enterobacteriaceae exclusivity organisms and
exclusive of the exclusivity strains commonly found
30 in stool due to the behavior of probe 1864. The expected behaviour of probe pairs of probe sets I

-39-

and II in a capture detection format with respect to EIEC, non-EIEC Enterobacteriaceae, and strains commonly found in stool is summarized in Table 13, in which a (+) signal is used as the cut-off for
05 detection. (In Table 13, as in Tables 11 and 12, any pair selected from probes 1500, 1501 and 1911 displays the behavior indicated.)

Assay Formats for Probes

Probes, probe pairs, and probe sets can be used
10 in a variety of hybridization assay formats. Such hybridization assays include solution hybridization assays in which the sequences to be detected and the probes are free in solution, or assays in which one of the sequence or probe is fixed to a solid
15 support. Shigella specific fragments, portions thereof, oligonucleotide probes derived from the fragments, complements, or homologues can be used in dot blot formats or other appropriate hybridization-based assay formats. For example, the
20 large fragments or portions thereof can be prepared as probes by nick translation or other suitable methods for filter hybridization (see e.g., U.S. Patent 4,358,535, Falkow et al.).

The probes can be used in suitable capture
25 detection assay formats (see e.g., D.V. Morrissey, et al., Analytical Biochemistry, 181:345-359 (1989); W.R. Hunsaker, et al., Analytical Biochemistry, 181:360-370 (1989); H. Lomeli, et al., Clinical Chemistry, 35:1826-1831 (1989); Pritchard, C.G. and
30 J.E. Stefano, Ann. Biol. clin. 48:492-497 (1990)).

-40-

In a capture/detector format, the probe pairs are preferably selected from non-overlapping portions of the same strand (a selected strand) of a Shigella specific fragment or a variant of a Shigella specific fragment. The probes can be separated by a distance consistent with activity in the selected assay. Thus, in a standard capture/detection format, the probes should be close enough that sample preparation does not separate the complementary sequences to the extent that the desired sensitivity of detection is compromised.

RNA probes may also be prepared. For example, probe nucleotide sequences can be incorporated into a fal-st MDV cDNA construct, and transcribed from linearized plasmid using T7 RNA polymerase. A detection probe prepared in this way can be used with one or more capture probes and amplified in Q β replicase system (Pritchard, C.G. and J.E. Stefano, Ann. Biol. Clin. 48: 492-497 (1990)).

The oligonucleotide probes described or others based on the sequence of the Shigella specific fragments can be used in the polymerase chain reaction. A second oligonucleotide can be prepared from the opposite strand.

The present invention will now be illustrated by the following examples, which are not intended to be limiting in any way.

-41-

EXAMPLESExample 1 Cyto-dot Panels

All cyto-dot panels were prepared by spotting approximately 1×10^8 cells of each bacterial isolate in a $5 \mu\text{l}$ volume onto nitrocellulose. The bacteria were lysed and the DNA was denatured by placing the nitrocellulose filters on 3MM paper wetted with 0.5 M NaOH and 1.5 M NaCl for 10 minutes. Following this treatment, the nitrocellulose filters were neutralized by placing them on 3MM paper wetted with 1 M Tris pH 7.5 and 1.5 M NaCl for 10 minutes. The latter neutralization step was repeated, and the DNA was fixed to the filters by baking under vacuum for 1-1.5 hours at 80°C .

Example 2 Hybridization Conditions

The hybridization conditions for all nick-translated fragments were as follows:
Prehybridization--was in 10X Denhardt's, 5X SET, 0.1 M phosphate buffer pH 7, 0.1% sodium pyrophosphate, 0.1% SDS for 3 hours at 65°C . (Note that 20X SET is 3M NaCl, 0.4M Tris-HCl pH 7.5 and 20mM EDTA).
Hybridization--was in 2X Denhardt's, 5X SET, 0.1 M phosphate buffer pH 7, 0.1% sodium pyrophosphate, 0.1% SDS and 1×10^6 counts of nick-translated probe per ml of hybridization solution. Hybridizations occurred overnight at 65°C . The autoradiographs were exposed for 15 hours.

The hybridization conditions for all kinased oligonucleotides (except probe 1864, a 17 base (b) oligonucleotide) were as follows:

-42-

Prehybridization--was in 5X Denhardt's, 6X SET, 0.1 M phosphate buffer pH 7, 0.1% sodium pyrophosphate, 0.1% SDS for 3 hours at 60 °C. Hybridization--was in 1X Denhardt's, 6X SET, 0.1 M phosphate buffer pH 7, 0.1% sodium pyrophosphate, 0.1% SDS and 1 X 10⁶ counts per minute of kinased oligonucleotide probe per ml of hybridization solution. Hybridizations occurred overnight at 60 °C. The autoradiographs were exposed for 15 hours or 7 days. The data recorded in Tables 2-8 are from 7 day exposures. The results of the two exposures were similar. The hybridization conditions for probe 1864 (17 b) were identical to those above except that the prehybridization, hybridization and wash temperatures were 50 °C rather than 60 °C.

-43-

TABLE 1

SIMPLIFIED OUTLINE OF SHIGELLA CLASSIFICATION

	SPECIES	SEROLOGICAL	SEROLOGICAL
		SUBGROUP	TYPE(S)
5	<i>S. dysenteriae</i>	A	1 through 10
	<i>S. flexneri</i>	B	1 through 6
	<i>S. boydii</i>	C	1 through 18
	<i>S. sonnei</i>	D	1

SUBSTITUTE SHEET

-44-

TABLE 2 : HYBRIDIZATION RESULTS FOR SHIGELLA DYSENTERIAE SEROTYPES

Genus species Strain ID, serotype	1911			NT			NT			NT			ompA			CR3		
	NT	1500	NT	1501	11-2	1682	1683	1708	1709	18-1a	1712	1713	19-2	1684	1707	1706	437	1864
<i>Shigella dysenteriae</i>																		
RF970 , 1	-	-	-	-	-	-	-	-	-	-	-	-	+/-	-	++++	++++	++++	+++
RF952 , 1	-	-	-	-	-	-	-	-	-	-	-	-	+/-	-	++++	++++	++++	+++
IG703 , 1	-	-	-	-	-	-	-	-	-	++++	++++	++++	+/-	-	++++	++++	++++	+++
IG704 , 1	-	-	-	-	-	-	-	-	-	++++	++++	++++	+/-	-	++++	++++	++++	+++
IG705 , 1	-	-	-	-	-	-	-	-	-	++++	++++	++++	+/-	-	++++	++++	++++	+++
IG710 , 1	-	-	-	-	-	-	-	-	-	++++	++++	++++	+/-	-	++++	++++	++++	+++
IG826 , 1	-	-	-	-	-	-	-	-	-	++++	++++	++++	+/-	-	++++	++++	++++	+++
IG828 , 1	-	-	-	-	-	-	-	-	-	++++	++++	++++	+/-	-	++++	++++	++++	-/+
IG774 , 2	-	-	++	++	++	++	++	++	++	++	++	++	+/-	-	++++	++++	++++	+++
IG725 , 2	++	+++	++	++	++	++	++	++	++	++	++	++	+/-	-	++++	++++	++++	+++
IG861 , 3	++++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	++++	+++
IG940 , 3	++++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	++++	+++
IG941 , 3	++++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	++++	+++
IG942 , 4	++++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	++++	+++
IG824 , 4	++++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	++++	+++
IG862 , 4	++++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	++++	+++
IG863 , 5	++++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	++++	+++

-45-

TABLE 2 : CONT'D

Genus species		1911		NT										NT		NT		ompA		CR3	
Strain ID,	serotype	NT	1500	6	1501	11-2	1682	1683	1708	1709	18-1a	1712	1713	NT	1684	1685	1707	1706	437	1864	
Shigella dysenteriae																					
IG864 ,	6	++++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	-	++++	++++	
IG865 ,	7	+++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	-	++++	++++	
IG866 ,	8a	+++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	-	++++	++++	
IG867 ,	9a	++++	++++	++++	+++	+++	+/	+-	++	-	-	-	-	+/-	-	-	-	-	++++	++++	
IG943 ,	9	++++	++++	+++	+++	+++	+/	+-	-	-	-	-	-	+/-	-	-	-	-	++++	++++	
IG944 ,	9	++++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	-	++++	++++	
IG868 ,	10	-	-	+++	+++	+++	+/	+-	++	-	-	-	-	+/-	-	-	-	-	++++	++++	
IG706		++++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	-	++++	++++	

SUBSTITUTE SHEET

Genus species		1911		NT					NT			NT			NT			CR3		
Strain ID,	serotype	NT	6	1500	11-2	1682	1683	1708	1709	1712	1713	1684	ompA	1707	1706	437	1864			
Shigella flexneri																				
IG879 ,	1a	-	-	-	-	-	-	-	-	++++	++++	+++	-	++++	-	+++	+++			
IG945 ,	1a	-	-	-	-	-	-	-	-	+++	++++	+++	-	++++	-	+++	+++			
IG946 ,	1b	-	-	-	-	-	-	-	-	++++	++++	+++	-	++++	-	+++	+++			
IG878 ,	1b	-	-	-	-	-	-	-	-	++++	+++	++	-	++++	-	+++	++			
IG818 ,	1b	-	-	-	-	-	-	-	-	++++	++++	+++	-	++++	-	+++	+++			
RF973 ,	2a	-	-	-	-	-	-	-	-	++++	++++	++++	-	++++	-	++++	++++			
IG716 ,	2b	+	-	-	-	-	-	-	-	++++	++++	++++	-	++++	-	++++	++++			
IG877 ,	3a	++	++	-	-	-	-	-	-	++++	++++	++++	-	++++	-	+++	+++			
IG947 ,	3a	-	-	-	+++	-	++	++	+++	++++	++++	+++	-	++++	-	+++	+++			
IG876 ,	3b	-	-	-	+++	-	++	++	+++	++++	++++	++++	-	++++	-	+++	+++			
IG875 ,	3c	-	-	-	+++	-	+	++	+++	++++	++++	+++	-	++++	-	+++	+++			
IG825 ,	3	+++	++	-	-	-	-	-	-	++++	++++	++++	-	++++	-	+++	+++			
IG874 ,	4a	-	-	-	-	-	-	-	-	++++	++++	++++	-	++++	-	+++	+++			
IG948 ,	4a	-	-	-	-	-	-	-	-	+++	++++	+++	-	++++	-	+++	+++			
IG873 ,	4b	-	-	-	-	-	-	-	-	++++	++++	++++	-	++++	-	+++	+++			
IG838 ,	4	-	-	-	-	-	-	-	-	++++	++++	+++	-	++++	-	+++	+++			
IG872 ,	5	-	-	-	-	-	-	-	-	-	-	+++	-	++++	-	+++	+++			

BNSDOCID: <WO 9303187A2>

-47-

TABLE 3 : CONT'D

Genus species Strain ID, serotype	NT		1911			NT			1708			1709			NT 18-1a			1713			NT 19-2			1684 1685			ompA 1707 1706			CR3 437 1864		
	6	NT	1500	1501	11-2	1682	1683	1708	1709	1712	1713	1712	1713	1712	1713	1712	1713	1712	1713	1712	1713	1712	1713	1712	1713	1712	1713	1712	1713	1712	1713	
Shigella flexneri																																
IG949 , 5	++	++	++	++	-	-	-	-	-	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	
IG823 , 6	++	++	+++	++	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
IG871 , 6	++++	++++	++++	++	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
IG950 , 6	+++	+++	+++	++	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
RF951	++	++	++	+++	+	+	+	+	+	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	
RF944	++	++	++	+++	+	+	+	+	+	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	
RF947	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
RF947	++	++	++	+++	+	+	+	+	+	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	
RF950	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
IG763	++	++	++	+++	+	+	+	+	+	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	
IG764	++	++	++	+++	+	+	+	+	+	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	
IG765	++	++	++	+++	+	+	+	+	+	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	
IG766	++	++	++	+++	+	+	+	+	+	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	
IG767	++	++	++	+++	+	+	+	+	+	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	
IG768	++	++	++	+++	+	+	+	+	+	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	
IG770	++	++	++	+++	+	+	+	+	+	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	
IG771	++	++	++	+++	+	+	+	+	+	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	

SUBSTITUTE SHEET

-48-

TABLE 3 : CONT'D

Genus species		1911		NT										NT 1712 1713				NT 1684				ompA				CR3			
Strain ID, serotype	NT 6	1501	1500	11-2	1682	1683	1708	1709	18-1a				19-2				1685	1707	1706	437 1864									
Shigella flexneri																													
IG772	++	++	++++	+	-		++	+++	++++	++++	++++	++++	++++	++++	++++	-	-	++++	+++	+++	+++	+++							
IG773	++	++	++++	+	-		++	+++	++++	++++	++++	++++	++++	++++	++++	-	-	++++	++++	++++	++++	+++							
IG775	-	-	-	-	-		-	-	++++	++++	++++	++++	++++	++++	++++	-	-	++++	+++	++++	++++	+++							
IG777	++	++	++++	+	-		++	+++	++++	++++	++++	++++	++++	++++	++++	-	-	++++	++++	++++	++++	+++							
IG778	++	++	++++	+	-		++	+++	++++	++++	++++	++++	++++	++++	++++	-	-	++++	++++	++++	++++	+++							
IG782	-	-	-	-	-		-	-	+++	++++	++++	++++	++++	++++	++++	-	-	++++	+++	++++	++++	+++							
IG724	-	-	++++	-	++		++	+++	++++	++++	++++	++++	++++	++++	++++	-	-	++++	+++	++++	++++	+++							
IG743	-	-	++++	-	-		++	+++	++++	++++	++++	++++	++++	++++	++++	-	-	++++	+++	++++	++++	+++							
IG744	-	-	-	-	-		-	-	++++	++++	++++	++++	++++	++++	++++	-	-	++++	+++	++++	++++	+++							
IG737	-	-	-	-	-		-	-	++++	++++	++++	++++	++++	++++	++++	-	-	++++	+++	++++	++++	+++							
IG735	-	-	-	-	-		-	-	++++	++++	++++	++++	++++	++++	++++	-	-	++++	+++	++++	++++	+++							
IG738	-	-	-	-	-		-	-	++++	++++	++++	++++	++++	++++	++++	-	-	++++	+++	++++	++++	+++							
IG736	-	-	+++	-	-		++	+++	++++	++++	++++	++++	++++	++++	++++	-	-	++++	+++	++++	++++	+++							
IG741	-	-	-	-	-		-	-	-	-	-	-	-	+++	+++	-	-	++++	+++	++++	++++	+++							
IG740	-	-	-	-	-		-	-	++++	++++	++++	++++	++++	++++	++++	-	-	++++	+++	++++	++++	+++							
IG739	-	-	-	-	-		-	-	++++	++++	++++	++++	++++	++++	++++	-	-	++++	+++	++++	++++	+++							
IG742	-	-	++++	-	++		++	+++	++++	++++	++++	++++	++++	++++	++++	-	-	++++	+++	++++	++++	+++							

-49-

TABLE 3 : CONT'D

Genus species Strain ID, serotype	1911			NT 1500 NT			1501 11-2			1682 1683 1708 1709			NT 1712 1713			NT 1684			ompA			CR3		
	6	NT	1500	11-2	1682	1683	1708	1709	18-1a	1712	1713	19-2	1685	1707	1706	437	1864							
<i>Shigella flexneri</i>																								
IG817	-	-	-	-	-	-	-	-	+++	++++	++++	+++	+++	-	-	+++	+++	+++	-	-	+++	+++	+++	+++
IG820	-	-	-	-	-	-	-	-	++++	++++	++++	+++	++	-	-	+++	+++	+++	-	-	+++	+++	+++	+++
IG822	-	-	-	-	-	-	-	-	+++	++++	++++	+++	+++	-	-	+++	+++	+++	-	-	+++	+++	+++	+++
IG709	-	-	-	-	-	-	-	-	-	-	-	+++	++++	-	-	+++	+++	+++	-	-	+++	+++	+++	+++
IG711	-	-	-	-	-	-	-	-	++++	-	++++	+++	+++	-	-	+++	+++	+++	-	-	+++	+++	+++	+++
IG726	++	++	+++	++	+/	++	++	++	++++	++++	++++	++++	++++	-	-	+++	+++	+++	-	-	+++	+++	+++	+++

SUBSTITUTE SHEET

-50-

TABLE 4 : HYBRIDIZATION RESULTS FOR SHIGELLA BOYDII SEROTYPES

Genus species Strain ID, serotype	1911			NT			NT			NT			ompA			CR3		
	NT	1500	1501	NT	11-2	1682	1683	1708	1709	NT	1712	1713	NT	1684	1707	1706	437	1864
	6	1501	11-2	1682	1683	1708	1709	18-1a					19-2	1685	1707	1706	437	1864
<i>Shigella boydii</i>																		
IG832 , 1	+++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	++++	++++
RF971 , 2	++++	++++	-	ND	ND	ND	ND	-	ND	ND	ND	ND	+/-	ND	ND	ND	++++	++++
IG880 , 3	++++	++++	-	-	-	-	-	-	-	-	-	-	-	-	-	-	++++	+++
IG935 , 4	++++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	++++	++++
IG882 , 5	-	-	+++	+++	++++	++	+++	-	++	-	-	-	++++	++++	++++	++++	++++	++++
IG936 , 5	-	-	+++	+++	++++	++	+++	-	++	-	-	-	++++	++++	++++	++++	++++	+++
IG883 , 6	++++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	++++	++++
IG884 , 7	+	-	+++	+++	++++	++	+++	-	++	-	-	-	++++	++++	++++	++++	++++	-/+
IG885 , 8	++++	++++	-	-	-	-	-	-	-	-	-	-	+	-	-	-	++++	++++
IG829 , 9	-	-	+++	+++	++++	++	+++	-	++	-	-	-	+++	+++	+++	+++	+++	+++
IG886 , 9	-	-	-	-	-	-	-	-	-	-	-	-	+++	+++	+++	+++	+++	+/-
IG887 , 10	++++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	++++	++++
IG937 , 10	++++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	++++	++++
IG3231, 10	++++	++++	-	-	-	-	-	-	-	-	-	-	+/-	-	-	-	++++	++++
IG938 , 11	-	-	+++	+++	++++	++	+++	-	++	-	-	-	+++	+++	+++	+++	+++	+++
IG888 , 11	-	-	+++	+++	++++	++	+++	-	++	-	-	-	+++	+++	+++	+++	+++	+/-
IG889 , 12	+	-	+++	+++	++++	-	+++	-	-	-	-	-	+/-	-	+++	+++	+++	+++

SUBSTITUTE SHEET

TABLE 4 : CONT'D

Genus species Strain ID, serotype	1911		NT		NT 1712 1713				NT 1684				ompA		CR3		
	NT 6	1500 1501	11-2	1501	18-1a	1709	1708	1683	1682	1683	1708	1709	1685	1707	1706	437	1864
Shigella boydii																	
RF974 , 13	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-/+
IG890 , 13	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
IG891 , 14	++++	++++	-	-	-	-	-	-	-	-	-	-	-	-	-	++++	++++
IG939 , 14	+++	++++	-	-	-	-	-	-	-	-	-	-	-	-	-	++++	++++
IG892 , 15	-	-	-	-	-	-	-	-	-	-	-	-	-	+++	+++	+++	+/-
IG700 , 16	-	-	+++	+++	-	++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++	+++
IG701 , 17	-	-	++++	+	-	+/-	+++	+++	+	+++	+++	+++	+++	++	+++	+++	+++
IG702 , 18	++++	++++	-	-	-	-	-	-	-	-	-	-	-	-	-	++++	++++
RF948	++++	++++	-	-	-	-	-	-	-	-	-	-	-	-	-	++++	++++
IG718	+++	++++	-	-	-	-	-	-	-	-	-	-	-	-	-	++++	+++

SUBSTITUTE SHEET

-52-

TABLE 5 : HYBRIDIZATION RESULTS FOR SHIGELLA SONNEI SEROTYPES

Genus species Strain ID	1911			NT			NT 1712 1713			NT 1684			ompA			CR3		
	NT 6	1500	1501	11-2	1682	1683	1708	1709	18-1a	1712	1713	NT 19-2	1685	1707	1706	437	1864	
<i>S. sonnei</i>																		
IG827	++++	++++	++++	+++	+++	++++	++	+++	-	-	-	+/-	-	-	-	+++	+++	+++
IG821	++++	++++	++++	+++	+++	++++	++	+++	-	-	-	+/-	-	-	-	+++	+++	+++
IG869	++++	++++	++++	+++	+++	+++	++	+++	-	-	-	+/-	-	-	-	+++	+++	+++
IG870	++++	++++	++++	+++	+++	+++	++	+++	-	-	-	+/-	-	-	-	+++	+++	+++
IG929	++++	++++	++++	+++	+++	++++	++	+++	-	-	-	+++	+++	-	-	+++	+++	+++
IG930	++++	++++	++++	+++	+++	++++	++	+++	-	-	-	+/-	-	-	-	+++	+++	+++
IG931	++++	++++	++++	+++	+++	++++	++	+++	-	-	-	+/-	-	-	-	+++	+++	+++
IG932	++++	++++	++++	+++	+++	++++	++	+++	-	-	-	+/-	-	-	-	+++	+++	+++
IG933	++++	++++	++++	+++	+++	++++	++	+++	-	-	-	+/-	-	-	-	+++	+++	+++
IG934	++++	++++	++++	+++	+++	++++	++	+++	-	-	-	+/-	-	-	-	+++	+++	+++
IG951	++++	++++	++++	+++	+++	++++	++	+++	-	-	-	+++	+++	-	-	+++	+++	+++
IG952	++++	++++	++++	+++	+++	++++	++	+++	-	-	-	+++	+++	-	-	+++	+++	+++
IG953	++++	++++	++++	+++	+++	++++	++	+++	-	-	-	+++	+++	-	-	+++	+++	+++
IG954	++++	++++	++++	+++	+++	++++	++	+++	-	-	-	+++	+++	-	-	+++	+++	+++
IG955	++++	++++	++++	+++	+++	++++	++	+++	-	-	-	+++	+++	-	-	+++	+++	+++
IG956	++++	++++	++++	+++	+++	++++	++	+++	-	-	-	+++	+++	-	-	+++	+++	+++
IG957	++++	++++	++++	+++	+++	++++	++	+++	-	-	-	+++	+++	-	-	+++	+++	+++

-53-

TABLE 5 : CONT'D

Genus species Strain ID	1911			NT			NT			NT			NT			ompA			CR3		
	NT	1500	1501	NT	11-2	1682	1683	1708	1709	NT	1712	1713	NT	1684	1685	1707	1706	437	1864	437	1864
<i>S. sonnei</i>																					
IG958	++++	++++	++++	+++	+++	+++	++++	++	+++	-	-	-	+++	+++	+++	-	-	+++	+++	+++	+++
IG959	++++	++++	++++	+++	+++	+++	++++	++	+++	-	-	-	+++	+++	+++	-	-	+++	+++	+++	+++
IG960	++++	++++	++++	+++	+++	+++	++++	++	+++	-	-	-	+++	+++	+++	-	-	+++	+++	+++	+++
IG961	++++	++++	++++	+++	+++	+++	++++	++	+++	-	-	-	+++	+++	+++	-	-	+++	+++	+++	+++
IG962	++++	++++	++++	+++	+++	+++	++++	++	+++	-	-	-	+++	+++	+++	-	-	+++	+++	+++	+++
IG963	++++	++++	++++	+++	+++	+++	++++	++	+++	-	-	-	+++	+++	+++	-	-	+++	+++	+++	+++
IG964	++++	++++	++++	+++	+++	++	++++	++	+++	-	-	-	+++	+++	+++	-	-	+++	+++	+++	+++
IG965	++++	++++	++++	+++	+++	+++	++++	+++	+++	-	-	-	+++	+++	+++	-	-	+++	+++	+++	+++
RF968	++++	++++	++++	+++	+++	+++	+++	++	+++	-	-	-	+++	+++	+++	-	-	+++	+++	+++	+++
RF943	++++	++++	++++	+++	+++	+++	++++	++	+++	-	-	-	+/	-	-	-	-	+++	+++	+++	+++
RF949	++++	++++	++++	+++	+++	+++	++++	++	+++	-	-	-	+/	-	-	-	-	+++	+++	+++	+++
IG781	++++	++++	++++	+++	+++	+++	++++	++	+++	-	-	-	+/	-	-	-	-	+++	+++	+++	+++
IG723	++++	++++	++++	+++	+++	+++	++++	+++	+++	-	-	-	+/	-	-	-	-	+++	+++	+++	+++
IG717	++++	++++	++++	+++	+++	+++	++++	++	+++	-	-	-	+/	-	-	-	-	+++	+++	+++	+++
IG713	++++	++++	++++	+++	+++	+	+++	++	+++	-	-	-	+/	-	-	-	-	+++	+++	+++	+++
IG714	++++	++++	++++	+++	+++	+++	++++	++	+++	-	-	-	+/	-	-	-	-	+++	+++	+++	+++
IG720	+++	+++	+++	+++	+++	+++	++++	++	+++	-	-	-	+/	-	-	-	-	+++	+++	+++	+++

SUBSTITUTE SHEET

TABLE 5 : CONT'D

Genus species		1911										NT					NT 1712 1713					NT 1684					ompA					CR3																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																													
Strain ID	NT 6	1500	1501	11-2	1682	1683	1708	1709	18-1a	1712	1713	19-2	1685	1707	1706	437	1864																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																												</

-55-

TABLE 5 : CONT'D

Genus species Strain ID	1911		NT		NT 1712 1713		NT 1712 1713		NT 1684		ompA		CR3	
	NT 6	1501	11-2	1682	1683	1708	1709	18-1a	19-2	1685	1707	1706	437	1864
<i>S. sonnei</i>														
IG967	++++	++++	+++	++	+++	++	+++	-	+++	+++	-	-	+++	+++
IG968	++++	++++	++++	+++	++++	+++	++++	-	+++	+++	-	-	+++	+++
IG969	++++	++++	+++	+++	++++	++	+++	-	+++	+++	-	-	+++	+++
IG970	++++	++++	++++	++++	++++	+++	++++	-	+++	+++	-	-	+++	+++
IG971	++++	++++	++++	++++	++++	+++	++++	-	+++	+++	-	-	+++	+++
IG972	++++	++++	+++	+++	++++	++	+++	-	+++	+++	-	-	+++	+++
IG974	++++	++++	++++	++++	++++	+++	++++	-	+++	+++	-	-	+++	+++
IG975	++++	++++	+++	++	+++	++	+++	-	+++	+++	-	-	+++	+++
IG976	++++	++++	++++	++++	++++	+++	++++	-	+++	+++	-	-	+++	+++
IG979	++++	++++	+++	+++	+++	++	+++	-	+++	+++	-	-	+++	+++
IG980A	++++	++++	+++	+++	++++	++	+++	-	+++	+++	-	-	+++	+++
IG980B	++++	++++	+++	+++	++++	++	+++	-	+++	+++	-	-	+++	+++
IG982	++++	++++	+++	+++	++++	++	+++	-	+++	+++	-	-	+++	+++

SUBSTITUTE SHEET

-56-

TABLE 6A : HYBRIDIZATION OF SHIGELLA PROBES TO E. COLI

Genus species	Strain ID	NT	1911	NT				
		6	1500 1501	11-2	1682	1683	1708	1709
Enteroinvasive								
E. coli								
"	3138	++++	++++	-	-	-	-	-
"	3145	++++	++++	+++	+	+/-	+++	++++
"	3146	++	++++	-	-	-	-	-
"	3157	++++	++++	-	-	-	-	-
"	3037	++++	++++	-	-	-	-	-
Enterotoxigenic								
E. coli								
"	3118	-	-	-	-	-	-	-
"	3119	-	-	-	-	-	-	-
"	3120	-	-	-	-	-	-	-
"	3123	-	-	-	-	-	-	-
"	3127	-	-	-	-	-	-	-
"	3129	-	-	-	-	-	-	-
"	3132	-	-	-	-	-	-	-
"	3134	-	-	-	-	-	-	-
"	3136	+	-	-	-	-	-	-
"	3142	-	-	-	-	-	-	-
"	3147	-	-	-	-	-	-	-
"	3151	-	-	-	-	-	-	-
"	3154	-	-	-	ND	-	-	-
"	3156	-	-	-	-	-	-	-
"	3158	+	-	-	-	-	-	-
"	3160	-	-	-	-	-	-	-

-57-

TABLE 6A : CONT'D

Genus species	Strain ID	NT 6	1911	NT				
			1500 1501	11-2	1682	1683	1708	1709
Enteropathogenic								
E. coli - O157 serotype								
"	3137	-	-	-	-	-	-	-
"	3140	-	-	-	-	-	-	-
"	3143	-	-	-	-	-	-	-
"	3144	-	-	-	-	-	-	-
"	3152	-	-	-	-	-	-	-
"	3155	-	-	-	-	-	-	-
"	3162	-	-	-	-	-	-	-
"	3163	-	-	-	-	-	-	-
"	3164	-	-	-	-	-	-	-
"	3165	-	-	-	-	-	-	-
"	3040	-	-	-	-	-	-	-
Enteropathogenic								
E. coli - non O157 serotypes								
"	3038	+	-	+++	+	+++	++	+++
"	3039	+	-	-	-	-	-	-
"	3041	+	-	-	-	-	-	-
"	3042	-	-	-	-	-	-	-
"	3043	-	-	-	-	-	-	-
"	839	-	-	-	-	-	-	-
"	840	+	-	-	-	-	-	-
"	841	-	-	-	-	-	-	-
"	842	-	-	-	-	-	-	-
"	843	-	-	+++	-	-	++	+++
"	844	-	-	-	-	-	-	-
"	845	-	-	-	-	-	-	-

-58-

TABLE 6A : CONT'D

Genus species	Strain ID	NT 6	1911	NT				
			1500 1501	11-2	1682	1683	1708	1709
Enteropathogenic								
E. coli non - O157 serotype								
"	846	++++	-	-	-	-	-	-
"	847	-	-	-	-	-	-	-
"	848	-	-	-	-	-	-	-
"	3149	-	-	-	-	-	-	-
Non-pathogenic								
E. coli								
"	3116	-	-	-	-	-	-	-
"	3117	-	-	-	-	-	-	-
"	3121	-	-	-	-	-	-	-
"	3122	-	-	-	-	-	-	-
"	3124	-	-	-	-	-	-	-
"	3125	-	-	-	-	-	-	-
"	3126	-	-	-	-	-	-	-
"	3128	-	-	-	-	-	-	-
"	3130	-	-	-	-	-	-	-
"	3131	-	-	-	-	-	-	-
"	3133	-	-	-	-	-	-	-
"	3135	-	-	-	-	-	-	-
"	3139	-	-	-	-	-	-	-
"	3148	-	-	-	-	-	-	-
"	3150	-	-	-	-	-	-	-
"	3153	-	-	-	-	-	-	-
"	3158	-	-	-	-	-	-	-
"	3161	-	-	-	-	-	-	-

-59-

TABLE 6B : HYBRIDIZATION OF SHIGELLA PROBES TO E. COLI

Genus species	Strain ID	NT 18-1a	1712	1713	NT 19-2	1684 1685	ompA 1707 1706	CR3 437 1864
Enteroinvasive E. coli								
"	3138	-	-	-	+/-	-	-	- +++++
"	3145	-	-	-	+/-	-	-	+ +++
"	3146	-	-	-	++++	++++	-	+ +++
"	3157	-	-	-	+/-	-	-	++ +++
"	3037	-	-	-	+/-	-	-	++ -/+
Enterotoxigenic E. coli								
"	3118	-	-	-	+/-	-	-	- -
"	3119	-	-	-	+/-	-	-	- -
"	3120	-	-	-	+/-	-	- -/+	- -/+
"	3123	-	-	-	+/-	-	- -/+	- -
"	3127	-	-	-	+/-	-	- -/+	- -
"	3129	-	-	-	+/-	-	- -/+	++ -/+
"	3132	-	-	-	+/-	-	-	+ -
"	3134	-	-	-	+/-	-	-	- -
"	3136	-	-	-	+/-	-	-	+ -
"	3142	-	-	-	+/-	-	-	- -
"	3147	-	-	-	+/-	-	- -/+	- -
"	3151	-	-	-	+/-	-	-	- -
"	3154	-	ND	ND	+/-	-	-	- -
"	3156	-	-	-	+/-	-	-	++ -/+
"	3158	-	-	-	+	-	-	++ ND
"	3160	-	-	-	+/-	-	-	+ -

-60-

TABLE 6B : CONT'D

Genus species	Strain ID	NT 18-1a	1712	1713	NT 19-2	1684 1685	ompA 1707 1706	CR3 437 1864
------------------	--------------	-------------	------	------	------------	--------------	-------------------	-----------------

Enteropathogenic

E. coli - O157 serotype

"	3137	-	-	-	+/-	-	-	-	-	-
"	3140	-	-	-	+/-	-	-	-	-	-
"	3143	-	-	-	+/-	-	-	-	-	-
"	3144	-	-	-	+/-	-	-	-	-	-
"	3152	-	-	-	+/-	-	-	-	-	-
"	3155	-	-	ND	+/-	-	-	-	-	-
"	3162	-	-	-	+/-	-	-	-	-	-
"	3163	-	-	-	+/-	-	-	-	-	-
"	3164	-	-	-	+/-	-	-	-	-	-
"	3165	-	-	-	+/-	-	-	-	-	-
"	3040	-	-	-	+/-	-	-	-	-	-

Enteropathogenic

E. coli - non O157 serotypes

"	3038	-	-	-	+/-	-	-	-	-	-
"	3039	+++	++++	++++	+/-	-	-	-	+	-
"	3041	-	-	-	+	-	-	-	++	ND
"	3042	-	-	-	+/-	-	-	-	-	-
"	3043	-	-	-	+/-	-	-	-	-	-
"	839	-	-	-	+/-	-	-	-	-	-
"	840	-	-	-	+/-	-	-	-	++	-
"	841	-	-	-	+/-	-	-	-	-	-
"	842	-	-	-	+/-	-	-	-/+	-	-
"	843	-	-	-	+/-	-	-	-	-	-
"	844	-	-	-	+/-	-	-	-	-	-/+
"	845	-	-	-	+/-	-	-	-	+	-

-61-

TABLE 6B : CONT'D

Genus species	Strain ID	NT 18-1a	1712	1713	NT 19-2	1684 1685	ompA 1707 1706	CR3 437 1864
Enteropathogenic E. coli non - O157 serotype								
"	846	-	-	-	+/-	-	-	-
"	847	-	-	-	+/-	-	-	-
"	848	-	-	-	+/-	-	-	-
"	3149	-	-	-	+/-	-	-	-
Non-pathogenic E. coli								
"	3116	-	-	-	+/-	-	-	++
"	3117	-	-	-	+/-	-	-	-
"	3121	-	-	-	+/-	-	-	-
"	3122	-	-	-	+/-	-	-	-
"	3124	-	-	-	+/-	-	-	-
"	3125	-	-	-	+/-	-	-	-
"	3126	-	-	-	+/-	-	-	-
"	3128	-	-	-	+/-	-	-	+
"	3130	-	-	-	+/-	-	-	-
"	3131	-	-	-	+/-	-	-/+	-
"	3133	-	-	-	+/-	-	-	-
"	3135	-	-	-	+/-	-	-	-
"	3139	-	-	-	+/-	-	-/+	-
"	3148	-	-	-	+/-	-	-	-
"	3150	-	-	-	+/-	-	-	+
"	3153	-	-	-	+/-	-	-/+	-
"	3158	-	ND	ND	+/-	-	-	+
"	3161	-	-	-	+/-	-	-/+	+

SUBSTITUTE SHEET

TABLE 7A : HYBRIDIZATION OF SHIGELLA PROBES TO MISC. ENTEROBACTERIACEAE

5	Genus species	GT#	ATCC#	1911				1708	1709
				NT	1500	NT	11-2		
				6	1501	11-2	1682	1683	
	Alteromonas putrefaciens	1495	8071	-	-	-	-	-	-
	Acinetobacter calcoaceticus	1972		-	-	-	ND	ND	ND
	Citrobacter diversus	1608		-	-	-	-	-	-
10	" diversus	1475	27156	-	-	-	-	-	-
	" amalonaticus	1607		-	-	-	-	-	-
	" amalonaticus	0689	25406	-	-	-	-	-	-
	" amalonaticus	0690	25405	-	-	-	-	-	-
	" freundii	0041		-	-	-	-	-	-
15	" freundii	0031		-	-	-	-	-	-
	" freundii	1597		-	-	-	-	-	-
	" freundii	1476	29935	-	-	-	-	-	-
	" freundii	1477	33128	-	-	-	-	-	-
	" freundii	1491	8090	-	-	-	-	-	-
20	" freundii	1591		-	-	-	-	-	-
	" freundii	1595		-	-	-	-	-	-
	" freundii	1599		-	-	-	-	-	-
	" freundii	0685		-	-	-	-	-	-
	" freundii	3241		-	-	-	-	-	-

TABLE 7A : CONT'D

	Genus species	GT#	ATCC#	1911		NT	1500	NT		11-2	1682	1683	1708	1709
				6	1501		1501	11-2	1682					
5	Citrobacter freundii	0038		-	-	-	-	-	-	-	-	-	-	-
	"	0036		-	-	-	-	-	ND	ND	ND	ND	ND	ND
	"	0035		-	-	-	-	-	-	-	-	-	-	-
	"	0034		-	-	-	-	-	-	-	-	-	-	-
	"	0033		-	-	-	-	-	-	-	-	-	-	-
10	"	0032		-	-	-	-	-	-	-	-	-	-	-
	"	0037		-	-	-	-	-	-	-	-	-	-	-
	"	0039		-	-	-	-	-	ND	-	-	-	-	-
	"	0040		-	-	-	-	-	-	-	-	-	-	-
	Enterobacter aerogenes	1487	29940	-	-	-	-	-	-	-	-	-	-	-
15	"	0047	13048	-	-	-	-	-	-	-	-	-	-	-
	"	0048		-	-	-	-	-	-	-	-	-	-	-
	"	0049		-	-	-	-	-	-	-	-	-	-	-
	"	1467	29917	-	-	-	-	-	-	-	-	-	-	-
	"	1468	29918	-	-	-	-	-	-	-	-	-	-	-
20	"	1469	29919	-	-	-	-	-	-	-	-	-	-	-
	"	1470	29920	-	-	-	-	-	-	-	-	-	-	-
	"	1471	29921	-	-	-	-	-	-	-	-	-	-	-
	"			-	-	-	-	-	-	-	-	-	-	-
	"			-	-	-	-	-	-	-	-	-	-	-

SUBSTITUTE SHEET

TABLE 7A : CONT'D

Genus species	GT#	ATCC#	1911		1911		1682	1683	1708	1709
			NT	1500	NT	1501				
5			6		11-2					
	Enterobacter agglomerans	1472	29922	-	-	-	-	-	-	-
	" agglomerans	1473	29923	-	-	-	-	-	-	-
	" agglomerans	1488	29904	-	-	-	-	-	-	-
10	" agglomerans	1489	29915	-	-	-	-	-	-	-
	" agglomerans	1490	29916	-	-	-	-	-	-	-
	" agglomerans	1474	27998	-	-	-	-	-	-	-
	" amnigenus	1482	33072	-	-	-	-	-	-	-
	" cloacae	0052		-	-	-	-	-	-	-
15	" cloacae	3042		-	-	-	-	-	-	-
	" cloacae	0050		-	-	-	-	-	-	-
	" cloacae	3041		-	-	-	-	-	-	-
	" cloacae	1159		-	-	-	-	-	-	-
	" cloacae	1337		-	-	-	-	-	-	-
20	" cloacae	1481	29941	-	-	-	-	-	-	-
	" cloacae	1492	13047	-	-	-	-	-	-	-
	" cloacae	3043		-	-	-	-	-	-	-
	" gergoviae	1486	33028	-	-	-	-	-	-	-
	" intermedium	0677	33110	-	-	-	-	-	-	-

TABLE 7A : CONT'D

	Genus species	GT#	ATCC#	1911			NT			NT		
				6	1501	1500	6	1501	11-2	1682	1683	1708 1709
5	Enterobacter sakazakii	0063		-	-	-	-	-	-	-	-	-
	" sakazakii	1483	29544	-	-	-	-	-	-	-	-	-
	" taylorae	1497	35317	-	-	-	-	-	-	-	-	-
	" taylorae	0065		-	-	-	-	-	-	-	-	-
	Escherichia blattae	1460		-	-	-	-	-	-	-	-	-
10	" fergusonii	1453		-	-	-	-	-	-	-	-	-
	" fergusonii	1459		-	-	-	-	-	-	-	-	-
	" hermannii	1216	33650	-	-	-	-	-	-	-	-	-
	" vulneri	1456		-	-	-	-	-	-	-	-	-
	" vulneri	1217	33821	-	-	-	-	-	-	-	-	-
15	Hafnia alvei	0241	29927	-	-	-	-	-	-	-	-	-
	" alvei	1153		-	-	-	-	-	-	-	-	-
	Klebsiella oxytoca	1606		-	-	-	-	-	-	-	-	-
	" oxytoca	1605		-	-	-	-	-	-	-	-	-
	" oxytoca	1503	13182	-	-	-	-	-	-	-	-	-
20	" ozaenae	1499	11296	-	-	-	-	-	-	-	-	-
	" planticola	1478	33531	-	-	-	-	-	-	-	-	-
	" pneumoniae	1150		-	-	-	-	-	-	-	-	-

SUBSTITUTE SHEET

TABLE 7A : CONT'D

	Genus species	GT#	ATCC#	1911		NT			
				NT	6	1500	11-2	1501	1709
5	Klebsiella pneumoniae	1500	13883	+	-	-	-	-	-
	" pneumoniae	1502	29939	-	-	-	-	-	-
	" pneumoniae	0252		-	-	-	-	-	-
10	" pneumoniae	1177		-	-	-	-	-	-
	" terrigena	1479	33257	-	-	-	-	-	-
	Morganella morganii	1147		-	-	-	-	-	-
	Proteus mirabilis	1148		-	-	-	ND	-	ND
	" mirabilis	1208		-	-	-	-	-	-
15	" mirabilis	1493	25933	-	-	-	ND	ND	ND
	" mirabilis	1496	29906	-	-	-	-	-	-
	" mirabilis	1501	7002	-	-	-	-	-	-
	" vulgaris	0370		-	-	-	-	-	-
	Providencia stuartii	3044		-	-	-	-	-	-
20	Pseudomonas aeruginosa	3045		-	-	-	-	-	-
	Salmonella typhimurium	0389	23566	-	-	-	ND	ND	ND
	Serratia marcescens	0392	29937	-	-	-	ND	-	-
	" marcescens	1151		-	-	-	-	-	-
	Yersinia enterocolitica	0424		-	-	-	-	-	-
25	" enterocolitica	3219		-	-	-	-	-	-

TABLE 7B : HYBRIDIZATION OF SHIGELLA PROBES TO MISC. ENTEROBACTERIACEAE

5	Genus species	GT#	ATCC#	NT 1712 1713		NT 1684		ompA		CR3	
				18-1a		19-2	1685	1707	1706	437	1364
10	Alteromonas putrefaciens	1495	8071	-	-	+/-	-	-	-	-	-
	Acinetobacter calcoaceticus	1972		-	ND	+/-	ND	ND	ND	-	-
	Citrobacter diversus	1608		-	-	+/-	-	-	-/+	-	-
	" diversus	1475	27156	-	-	+/-	-	-	-/+	-	-
	" amalonaticus	1607		-	-	+/-	-	-	-/+	-	-
	" amalonaticus	0689	25406	-	-	+/-	-	-	-/+	-	-
	" amalonaticus	0690	25405	-	-	+/-	-	-	-/+	-	-
	" freundii	0041		-	-	+/-	-	-	-	-	-
	" freundii	0031		-	-	+/-	-	-	-	-	-
	" freundii	1597		-	ND	+/-	-	-	-	-	-
15	" freundii	1476	29935	-	-	+/-	-	-	-	-	-
	" freundii	1477	33128	-	-	+/-	-	-	-	-	-
	" freundii	1491	8090	-	-	+/-	-	-	-	-	-
	" freundii	1591		-	-	+/-	-	-	-	-	-
	" freundii	1595		-	-	+/-	-	-	-	-	-
20	" freundii	1599		-	-	+/-	-	-	-	-	-
	" freundii	0685		-	-	+/-	-	-	-	-	-
	" freundii	3241		-	-	+/-	-	-	-/+	-	-
	" freundii			-	-	+/-	-	-	-	-	-
	" freundii			-	-	+/-	-	-	-	-	-

SUBSTITUTE SHEET

SUBSTITUTE SHEET

TABLE 7B : CONT'D

Genus species	GT#	ATCC#	NT 1712 1713		NT 1684 19-2 1685	OmpA		CR3
			18-1a			1707	1706	
5								
Enterobacter	1472	29922	-	-	+/-	-	-	-
"	1473	29923	-	-	+/-	-	-	-
"	1488	29904	-	-	+/-	-	-	-
"	1489	29915	-	-	+/-	-	+/-	-
10	1490	29916	-	-	+/-	-	-	-
"	1474	27998	-	-	+/-	-	-	-
"	1482	33072	-	-	+/-	-	-	-
"	0052		-	-	+/-	-	+	-
"	3042		-	-	+/-	-	-	-
15	0050		-	-	+/-	-	-	-
"	3041		-	-	+/-	-	-	-
"	1159		-	-	+/-	-	-	-
"	1337		-	-	+/-	-	-	-
"	1481	29941	-	-	+/-	-	-	-
20	1492	13047	-	-	+/-	-	-	-
"	3043		-	-	+/-	-	-	-
"	1486	33028	-	-	+/-	-	-	-
"	0677	33110	-	-	+/-	-	-	-

TABLE 7B : CONT'D

Genus species	GT#	ATCC#	NT 1712 1713 18-1a	NT 1684 19-2 1685	ompA 1707 1706	CR3 437 1864
5						
Enterobacter sakazakii	0063		-	-	-	-
" sakazakii	1483	29544	-	-	-	-
" taylora	1497	35317	-	-	-	-
" taylora	0065		-	-	-	-
10 Escherichia blattae	1460		-	-	++++	-
" fergusonii	1453		-	-	+	-
" fergusonii	1459		-	-	++	-
" hermanii	1216	33650	-	ND	-	-
" vulneri	1456		-	-	-	-
15 " vulneri	1217	33821	-	-	-	-
Hafnia alvei	0241	29927	-	-	-	-
" alvei	1153		-	-	-	-
Klebsiella oxytoca	1606		-	-	-	ND
" oxytoca	1605		-	ND	-	ND
20 " oxytoca	1503	13182	-	-	-	-
" ozaenae	1499	11296	-	-	-	-
" planticola	1478	33531	-	-	-	-
" pneumoniae	1150		-	-	-	-

Genus species	GT#	ATCC#	NT 1712 1713 18-1a	NT 1684 19-2 1685	OmpA 1707 1706	CR3 437 1864
Klebsiella pneumoniae	1500	13883	-	-	-	-
" pneumoniae	1502	29939	-	+/-	-	-
" pneumoniae	0252		-	+/-	-	-
" pneumoniae	1177		-	+/-	-	-
" terrigena	1479	33257	-	+/-	-	-
Morganella morganii	1147		-	+/-	-	-
Proteus mirabilis	1148		-	+/-	-	-
" mirabilis	1208		-	+/-	-	-
" mirabilis	1493	25933	-	+/-	ND	-
" mirabilis	1496	29906	-	+/-	-	-
" mirabilis	1501	7002	-	+/-	-	-
" vulgaris	0370		-	+/-	-	-
Providencia stuartii	3044		-	+/-	-	-
Psuedomonas aeruginosa	3045		-	+/-	-	-/+
Salmonella typhimurium	0389	23566	-	+/-	ND	-
Serratia marcescens	0392	29937	-	+/-	-	-
" marcescens	1151		-	+/-	-	-
Yersinia enterocolitica	0424		-	+/-	-	-
" enterocolitica	3219		-	+/-	-	-

-72-

TABLE 8 : HYBRIDIZATION OF SHIGELLA PROBES TO BACTERIA COMMONLY FOUND IN STOOL

	Genus/Species	GT#	ATCC#	NT 6			NT 19-2			ompA			CR3
				1500	1501	1911	1684	1685	1706	1707	1864		
5	Acinetobacter calcoaceticus	0002	19606	-	-/+	-	-	-	-	-	-	-	
	Acinetobacter lwoffii	0004	9957	-	-	-	-	-	-	-	-	-	
	Aeromonas hydrophila	0006	7965	-	-/+	-	-	-	-	-	-	-	
	Aeromonas sobria	0007	9071	-	-	-	-	-	-	-	-	-	
	Alteromonas putrefaciens	1495	8071	-	-	-	-	-	-	-	-	-	
10	Citrobacter amalonaticus	0690	25405	-	-	-	-	-	-	-	-	-	
	Edwardsiella tarda	0569	15947	-	-	-	-	-	-	-	-	-	
	Haemophilus influenzae	0244	19418	-	+/-	-	-	-	+/-	-	-	-	
	Plesiomonas shigelloides	2197	14029	-	-	-	-	-	-	-	-	-	
	Providencia alcalifaciens	0371	9886	-	-/+	-	-	-	-	-	-	-	
15	Providencia rettgeri	0373	9944	-	-	-	-	-	-	-	-	-	
	Providencia stuartii	0375	29914	-	-	-	-	-	-	-	-	-	
	Salmonella arizona	0799	13314	-	-	-	-	-	-	-	-	-	
	Salmonella typhimurium	0389	23566	-	+/-	-	-	-	-	-	-	-	
	Vibrio parahemolyticus	0568	17802	-	-	-	-	-	-	-	-	-	
20	Xanthomonas maltophila	0417	13637	-	-	-	-	-	-	-	-	-	
	Yersinia enterocolitica	0419	9610	-	-	-	-	-	-	-	-	-	
	Yersinia pseudotuberculosis	0519	29833	-	-	-	-	-	-	-	-	-	

TABLE 8 : CONT'D

	Genus/Species	GT#	ATCC#	NT 6			NT 19-2	ompA		CR3
				1500	1501	1911		1706	1707	
5	Versinia ruckeri	0522	29473	-	+/-	-	-	-	-	-
	Alcaligenes denitrificans ss xyloxidans	0001	27062	-	-/+	-	-	-	-	-
	Alcaligenes faecalis	0610	875	-	-	-	-	-	-	-
10	Kingella kingae	0247	23330	-	-	-	-	-	-	-
	Kingella denitrificans	0245	33394	-	-	-	-	-	-	-
	Kingella indologenes	0246	25869	-	-	-	-	-	-	-
	Moraxella osloensis	0301	19962	-	-	-	-	-	-	-
	Neisseria cinerea	0307	14685	-	-	-	-	-	-	-
15	Neisseria flavescens	0310	13120	-	-/+	-	-	-	-	-
	Neisseria gonorrhoeae, type	0315	19424	-	-/+	-	-	-	-	-
	Neisseria meningitidis	0348		-	-/+	-	-	-	-	-
	Neisseria mucosa	0353	19696	-	-	-	-	-	-	-
	Pseudomonas acidovorans, type	0376	15668	-	-	-	-	-	-	-
20	Achromobacter xerosis	0810	14780	-	-	-	-	-	-	-
	Gardnerella vaginalis	0240	1408	-	-	-	-	-	-	-
	Acidaminococcus fermentans	3129	25085	-	-	-	-	-	-	-
	Ruminococcus bromeli	3142	27255	-	-	-	-	-	-	-

-74-

TABLE 8 : CONT'D

	Genus/Species	GT#	ATCC#	NT 6			NT 19-2	ompA		CR3
				1500	1501	1911		1706	1707	
5	Bacteroides gracilis	0716	33236	-	-	-	-	-	-	-
	Bacteroides ureolyticus	0715	33387	-	-/+	-	-	-	-	-
	Campylobacter jejuni, type	0022	33560	-	-	-	-	-	-	-
	Campylobacter coli	0016	33559	-	-/+	-	-	-	-	-
	Campylobacter laridis	0024	35223	-	-	-	-	-	-	-
10	Wolinella curva	2224	35224	-	-	-	-	-	-	-
	Wolinella recta	0718	33238	-	-	-	-	-	-	-
	Wolinella succinogenes	0614	29543	-	-	-	-	-	-	-
	Bacillus cereus	0008	14579	-	-	-	-	-	-	-
	Butyrivibrio fibrosolvens	3139	19171	-	-	-	-	-	-	-
15	Clostridium difficile	0043	9689	-	-/+	-	-	-	-	-
	Clostridium perfringens	0044	3624	-	-/+	-	-	-	-	-
	Clostridium sordellii	0567	9714	-	-/+	-	-	-	-	-
	Eubacterium lentum	2196	25559	-	-	-	-	-	-	-
	Eubacterium rectale	0236	35183	-	-/+	-	-	-	-	-
20	Lactobacillus acidophilus	0256	4356	-	-	-	-	-	-	-
	Lactobacillus casei	0805	393	-	-	-	-	-	-	-
	Lactobacillus minutus	0257	33267	-	-/+	-	-	-	-	-

-75-

TABLE 8 : CONT'D

Genus/Species	GT#	ATCC#	NT 6			NT 19-2			ompA		CR3
			1500	1501	1911	1684	1685	1706	1707	1864	
5											
Lactobacillus plantarum	0258	8014	-	-	-	-	-	-	-	-	-
Listeria grayi	0674		-	+/-	-	-	-	-	-	-	-
Listeria innocua	0260		-	+/-	-	-	-	-	-	-	-
Listeria innocua	0750		-	-/+	-	-	-	-	-	-	-
10											
Listeria ivanovii	1037		-	+/-	-	-	-	-	-	-	-
Listeria monocytogenes	1016		-	+/-	-	-	-	-	-	-	-
Listeria seeligeri	0287		-	-	-	-	-	-	-	-	-
Listeria welshrimpi	0291		-	+/-	-	-	-	-	-	-	-
Peptococcus asaccharolyticus	0360	29743	-	-	-	-	-	-	-	-	-
15											
Peptococcus magnus	0361	29328	-	-	-	-	-	-	-	-	-
Peptostreptococcus anaerobius	0359	27337	-	-	-	-	-	-	-	-	-
Sarcina maxima	0391	33910	-	-/+	-	-	-	-	-	-	-
Staphylococcus aureus	0399	12600	-	-/+	-	-	-	-	-	-	-
Staphylococcus epidermidis	0401	14990	-	-/+	-	-	-	-	-	-	-
20											
Streptococcus agalactiae	0405	13813	-	-/+	-	-	-	-	-	-	-
Streptococcus faecalis	0406	19433	-	-/+	-	-	-	-	-	-	-
Streptococcus faecium	0407	6056	-	-/+	-	-	-	-	-	-	-
Streptococcus mutans	0412	25175	-	-/+	-	-	-	-	-	-	-
Streptococcus salivarius	0410	13419	-	-/+	-	-	-	-	-	-	-

SUBSTITUTE SHEET

-76-

TABLE 8 : CONT'D

	Genus/Species	GT#	ATCC#	NT 6			NT 19-2			ompA	CR3
				1500	1501	1911	1684	1685	1706		
5	<i>Streptococcus sanguis</i>	0411	10556	-	+/-	-	-	-	-	-	-
	<i>Actinomyces israelii</i>	0005	10049	-	-	-	-	-	-	-	-
	<i>Bifidobacterium dentium</i>	0012	27534	-	-	-	-	-	-	-	-
10	<i>Bifidobacterium bifidum</i>	0571	35914	-	-	-	-	-	-	-	-
	<i>Corynebacterium genitalium</i>	0045	33031	-	-	-	-	-	-	-	-
	<i>Mycobacterium smegmatis</i>	0306	14468	-	-	-	-	-	-	-	-
	<i>Propionibacterium acnes</i>	0363	6919	-	-	-	-	-	-	-	-
	<i>Fusobacterium mortiferum</i>	0573	9817	-	-/+	-	-	-	-	-	-
15	<i>Fusobacterium necrophorum</i>	0238	25286	-	-	-	-	-	-	-	-
	<i>Mobiluncus mulieris</i>	0300	35243	-	-	-	-	-	-	-	-
	<i>Veionella atypica</i>	0413	14894	-	-	-	-	-	-	-	-
	<i>Bacteroides fragilis</i>	0010	23745	-	-/+	-	-	-	-	-	-
	<i>Bacteroides thetaiotaomicron</i>	0572	29741	-	-/+	-	-	-	-	-	-
20	<i>Bacteroides melaninogenicus</i>	0011	25845	-	-	-	-	-	-	-	-
	<i>Flavobacterium meningosepticum</i>	0237	13253	-	-/+	-	-	-	-	-	-
	<i>Candida albicans</i>	0028	18804	-	-	-	-	-	-	+/-	-
	<i>Candida glabrata</i>	0029	2001	-	-	-	-	-	-	-	-
	<i>Candida stellatoidae</i>	0609	36232	-	-	-	-	-	-	-	-
25	<i>Candida tropicalis</i>	0570	750	-	-	-	-	-	-	-	-

TABLE 9
LIST OF OLIGONUCLEOTIDE PROBES DESIGNED FROM SHIGELLA SPECIFIC FRAGMENTS

Parental DNA Fragment	Oligonucleotide Probe No. (SEQ ID NO)	Length of Probes (bases)	Sequence of oligonucleotide probes
NT 6 and adjacent 3' sequence (164 b)	1500 (SEQ ID NO:14)	35	5' TTGCAGCGCCTCTACTACGGGATACAGCCTCCATT 3'
	1501 (SEQ ID NO:15)	35	5' CCTCCTTCAGGGCGGATTCAGCCGTTACATTGT 3'
	1911 (SEQ ID NO:16)	40	5' CCGATCTTCTATTGTACGACGGTTCGTCAAAAGCTAAT 3'
NT 11-2 (796 b)	1682 (SEQ ID NO:17)	41	5' CTGGTGAACACGCTCTTACAAAGATGGTTCCTGGATGGATT 3'
	1683 (SEQ ID NO:18)	41	5' AGTCTTTCCGTGTTTCTCAGAAATGGGGCAACGTGCAAAA 3'
	1708 (SEQ ID NO:19)	35	5' CCACCGTTGAAGCGTAAACCGTTTGACCGATGGAT 3'
	1709 (SEQ ID NO:20)	36	5' GCTGGGGTCTACAGGTGCAATAACCACTTAGACGGT 3'
NT 14 (786 b)	437 (SEQ ID NO:21)	48	5' CGATGATGCCATTCTCTGCCAGCTCCGTCTGGAGCCGCCG- GGTTCC 3'
NT 15 (587 b)	1864 (SEQ ID NO:22)	17	5' GGAGCAGTCTGGTCTGA 3'
NT18-1a (630 b)	1712 (SEQ ID NO:23)	37	5' CCTGTGGCTCTCGGTTCTGATGGTATAGCAACTAAAT 3'
	1713 (SEQ ID NO:24)	37	5' CAAGGATTTTCGGGAATTGAGTGGGGAGTTGCGAAAT 3'
NT19-2 (388 b)	1684 (SEQ ID NO:25)	35	5' CAGGCAATCGAAGCATATCGCGGTTCTCCACAACT 3'
	1685 (SEQ ID NO:26)	35	5' TGAATGCGCTGACCGAAAAACAGCGCTGGGTATCT 3'
ompA	1706 (SEQ ID NO:27)	35	5' GTGATGGCCCCATTCAACACCACCTGCGAATACCGG 3'
	1707 (SEQ ID NO:28)	35	5' CTCAGATTACCTGTACATTGTTGTGAGCTTTGG 3'

-78-

Table 10
Summary of Probe Hybridization Results to Shigella Serotypes
(+ signal used as a cut-off for detectability)

SEROTYPE	Tested	NT-6	1911 1500 1501	NT11-2	1682	1683	1708	1709	NT18-1a	1712	1713	NT19-2	1684	1685	1707	1706	437	1864
<i>S. Dysenteriae</i>																		
1	8	-	-	-	-	-	-	-	6	6	6	-	-	-	8	8	7	
2	2	1	1	2	2	2	-	-	-	-	-	-	-	-	2	2	2	
3	3	3	3	-	-	-	-	-	-	-	-	-	-	-	-	-	3	3
4	3	3	3	-	-	-	-	-	-	-	-	-	-	-	-	-	3	3
5	1	1	1	-	-	-	-	-	-	-	-	-	-	-	-	-	1	1
6	1	1	1	-	-	-	-	-	-	-	-	-	-	-	-	-	1	1
7	1	1	1	-	-	-	-	-	-	-	-	-	-	-	-	-	1	1
8	1	1	1	-	-	-	-	-	-	-	-	-	-	-	-	-	1	1
9	3	3	3	2	2	-	1	-	-	-	-	-	-	-	-	-	3	3
10	1	-	-	1	1	-	1	-	-	-	-	-	-	-	-	-	1	1
untyped	1	1	1	-	-	-	-	-	-	-	-	-	-	-	-	-	1	1
	25																	
<i>S. flexneri</i>																		
1	5	-	-	-	-	-	-	-	5	5	5	5	5	5	-	-	5	5
2	2	1	-	-	-	-	-	-	2	2	2	2	2	2	-	-	2	2
3	5	2	2	3	-	3	3	3	5	5	5	5	5	5	-	-	5	5
4	4	-	-	-	-	-	-	-	4	4	4	4	4	4	-	-	4	4
5	2	1	1	-	-	-	-	-	1	1	1	2	2	2	-	-	2	2
6	3	3	3	-	-	-	-	-	-	-	-	-	-	-	-	-	3	3
untyped	36	16	16	20	16	2	20	20	34	33	34	36	36	36	-	-	36	36
	57																	

-79-

Table 10 (cont')

SEROTYPE	Tested	NT-6	1911	NT11-2	1682	1683	1708	1709	NT18-1a	1712	1713	NT19-2	1684	1685	1707	1706	437	1864
			1500															
			1501															
<hr/>																		
<i>S. boydii</i>																		
1	1	1	1	-	-	-	-	-	-	-	-	-	-	-	-	-	1	1
2	1	1	1	-	-	-	-	-	-	-	-	-	-	-	-	-	1	1
3	1	1	1	-	-	-	-	-	-	-	-	-	-	-	-	-	1	1
4	1	1	1	-	-	-	-	-	-	-	-	-	-	-	-	-	1	1
5	2	-	-	2	2	2	2	2	-	-	-	2	2	2	2	2	2	2
6	1	1	1	-	-	-	-	-	-	-	-	-	-	-	-	-	1	1
7	1	1	1	1	1	1	1	1	-	-	-	1	1	1	1	1	1	1
8	1	1	1	-	-	-	-	-	-	-	-	1	-	-	-	-	1	1
9	2	-	-	1	1	1	1	-	-	-	-	2	2	2	2	2	2	1
10	3	3	3	-	-	-	-	-	-	-	-	-	-	-	-	-	3	3
11	2	-	-	2	2	2	2	2	-	-	-	2	2	2	2	2	2	1
12	1	1	-	1	1	1	-	-	-	-	-	-	-	-	1	1	1	1
13	2	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	2	2	2	-	-	-	-	-	-	-	-	-	-	-	-	-	2	2
15	1	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	1	1
16	1	-	-	1	1	1	1	1	-	-	-	1	1	1	1	1	1	1
17	1	-	-	1	1	1	-	1	-	-	-	1	1	1	1	1	1	1
18	1	1	1	-	-	-	-	-	-	-	-	-	-	-	-	-	1	1
untyped	2	2	2	-	-	-	-	-	-	-	-	-	-	-	-	-	2	2
27																		
<hr/>																		
<i>S. sonnei</i>																		
1	64	64	64	64	64	64	64	64	-	-	-	33	33	33	-	-	64	64

-80-

TABLE 11 : SUMMARY OF HYBRIDIZATION OF SHIGELLA SPECIFIC FRAGMENTS AND
CLIGONUCLEOTIDE CAPTURE/DETECTION PROBES TO SHIGELLA SEROTYPES
(++ SIGNAL USED AS CUT-OFF FOR DETECTABILITY).

5	SEROTYPE	# TESTED	1911		NT11		1682		1708		NT18		1712		NT19		1684		ompA		CR3		
			NT6	1500	NT11	-2	1683	1709	-1a	1713	-2	1685	1707	1864									
10	S. dysenteriae																						
	1	8	-	-	-	-	-	-	6	6	-	-	-	8	7								
	2	2	1	1	2	2	-	-	-	-	-	-	2	2									
	3	3	3	3	-	-	-	-	-	-	-	-	-	3									
	4	3	3	3	-	-	-	-	-	-	-	-	-	3									
15	5	1	1	1	-	-	-	-	-	-	-	-	-	1									
	6	1	1	1	-	-	-	-	-	-	-	-	-	1									
	7	1	1	1	-	-	-	-	-	-	-	-	-	1									
	8	1	1	1	-	-	-	-	-	-	-	-	-	1									
	9	3	3	3	2	-	-	-	-	-	-	-	-	3									
20	10	1	-	-	1	-	-	-	-	-	-	-	-	1									
	untyped	1	1	1	-	-	-	-	-	-	-	-	-	1									
			25																				

-81-

TABLE 11 : CONT'D

5	SEROTYPE	# TESTED	1911										ompA		CR3
			NT6	1500	NT11	1682	1708	NT18	1712	NT19	1684	1706	437C		
				1501	-2	1683	1709	-1a	1713	-2	1685	1707	1864		
S. flexneri															
	1	5	-	-	-	-	-	5	5	5	5	-	-	5	
	2	2	-	-	-	-	-	2	2	2	2	-	-	2	
10	3	5	2	2	3	-	3	5	5	5	5	-	-	5	
	4	4	-	-	-	-	-	4	4	4	4	-	-	4	
	5	2	1	1	-	-	-	1	1	2	2	-	-	2	
	6	3	3	3	-	-	-	-	-	-	-	-	-	3	
	untyped	36	16	16	20	-	20	34	33	36	36	-	-	36	
15		57													
S. boydii															
	1	1	1	1	-	-	-	-	-	-	-	-	-	1	
	2	1	1	1	-	-	ND	-	ND	-	ND	ND	ND	1	
20	3	1	1	1	-	-	-	-	-	-	-	-	-	1	
	4	1	1	1	-	-	-	-	-	-	-	-	-	1	
	5	2	-	-	2	2	2	-	-	2	2	2	2	2	
	6	1	1	1	-	-	-	-	-	-	-	-	-	1	
	7	1	-	-	1	1	1	-	-	1	1	1	1	-	

-82-

TABLE 11 : CONT'D

5	SEROTYPE	# TESTED	1911		NT6	1500	1501	NT11		1682		1708		NT18		1712		NT19		ompA		CR3	
			-2	1683				1709	-1a	1713	-2	1685	1707	1864									
S. boydii																							
10	8	1	1	1	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	1	
	9	2	-	-	1	1	-	-	-	-	-	-	-	2	2	2	2	2	2	2	2	1	
	10	3	3	3	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	3	
	11	2	-	-	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	2	1	
	12	1	-	-	1	1	1	-	-	-	-	-	-	-	-	-	-	-	-	1	1	1	
	13	2	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
	14	2	2	2	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	2	
	15	1	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	1	1	-	
	16	1	-	-	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	
	17	1	-	-	1	1	-	-	-	-	-	-	-	-	1	1	1	1	1	1	1	1	
20	18	1	1	1	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	1	
	untyped	2	2	2	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	2	
		27																					
S. sonnei																							
	1	64	64	64	64	64	64	64	64	64	64	64	64	-	-	-	-	33	33	-	-	64	

TABLE 12

Summary of Inclusivity Hybridization Results of Shigella Specific Fragments and Capture/Detection Oligonucleotide Probes Included in Probe Sets 1 and 1 to Shigella Serotypes (++) SIGNAL USED AS CUT-OFF FOR DETECTABILITY).

- 5 Probe Set 1 includes oligonucleotides 1911/1500/1501, 1684/1685 and 1706/1707.
 Probe Set 2 includes oligonucleotides 1706/1707 and 1864/437c.

10	SEROTYPE	TESTED	NT6		1911	NT19		ompA		CR3		probe set I	probe set II	
					1500	1684	-2	1706	1707	437c				
					1501	1685			1864					
S. dysenteriae														
15	1	8	-	-	-	-	-	8	7	-	8	8		
	2	2	1	1	-	-	-	2	2	-	2	2		
	3	3	3	3	-	-	-	-	3	-	3	3		
	4	3	3	3	-	-	-	-	3	-	3	3		
	5	1	1	1	-	-	-	-	1	-	1	1		
20	6	1	1	1	-	-	-	-	1	-	1	1		
	7	1	1	1	-	-	-	-	1	-	1	1		
	8	1	1	1	-	-	-	-	1	-	1	1		
	9	3	3	3	-	-	-	-	3	-	3	3		
	10	1	-	-	-	-	-	-	1	-	-	1		
untyped		1	1	1	-	-	-	-	1	-	1	1		
TOTAL		25											24	25

TABLE 12 : CONT'D

5	SEROTYPE	TESTED	1911		OMPA			CR3		probe set I	probe set II
			NT6	1500 1501	NT19 -2	1684 1685	1706 1707	437c 1864			
S. flexneri											
10	1	5	-	-	5	5	-	5	5	5	
	2	2	-	-	2	2	-	2	2	2	
	3	5	2	2	5	5	-	5	5	5	
	4	4	-	-	4	4	-	4	4	4	
	5	2	1	1	2	2	-	2	2	2	
	6	3	3	3	-	-	-	3	3	3	
15	untyped	36	16	16	36	36	-	36	36	36	
	TOTAL	57							57	57	
S. boydii											
20	1	1	1	1	-	-	-	1	1	1	
	2	1	1	1	-	ND	ND	1	1	1	
	3	1	1	1	-	-	-	1	1	1	
	4	1	1	1	-	-	-	1	1	1	
	5	2	-	-	2	2	2	2	2	2	
	6	1	1	1	-	-	-	1	1	1	
	7	1	-	-	1	1	1	-	1	1	

TABLE 12 : CONT'D

5	SEROTYPE	TESTED	1911			ompA			CR3			probe set I	probe set II
			NT6	1500	1501	NT19 -2	1684	1706	1707	437C	1864		
S. boydii													
10	8	1	1	1	-	-	-	1	1	1	1	1	1
	9	2	-	-	2	2	2	1	1	1	2	2	2
	10	3	3	3	-	-	-	3	3	3	3	3	3
	11	2	-	-	2	2	2	1	1	1	2	2	2
	12	1	-	-	-	-	1	1	1	1	1	1	1
	13	2	-	-	-	-	-	-	-	-	-	-	-
	14	2	2	2	-	-	-	2	2	2	2	2	2
	15	1	-	-	-	-	1	-	-	1	1	1	1
	16	1	-	-	1	1	1	1	1	1	1	1	1
	17	1	-	-	1	1	1	1	1	1	1	1	1
20	18	1	1	1	-	-	-	1	1	1	1	1	1
	untyped	2	2	2	-	-	-	-	-	2	2	2	2
	TOTAL	27								25	25	25	25
S. sonnei													
	1	64	64	64	33	33	-	64	64	64	64	64	64

SUBSTITUTE SHEET

TABLE 13

SUMMARY OF HYBRIDIZATION RESULTS OF SHIGELLA PROBES, WHICH ARE INCLUDED IN PROBE SETS 1 AND 2, TO ENTEROINVASIVE E. COLI AND COMPETITOR ORGANISMS FOUND IN STOOL (+ SIGNAL USED AS CUT-OFF FOR DETECTABILITY).

	FRAGMENT (OLIGONUCLEOTIDES)	E. COLI (from Table 6)			OTHER ORGANISMS FOUND IN STOOL (from Tables 7,8)	
		EIEC 5	EVEC 43	Non- EVEC 18	Cyto-dot 91	DNA-dot 91
5	PROBE SET 1					
10	NT6 (1500/1501/1911)	5	-	-	-	-
15	NT19-2 (1684/1685)	1	-	-	-	-
	OMP A (1706/1707)	-	-	-	2	-
	PROBE SET 2					
	OMP A (1706/1707)	-	-	-	2	-
20	Class 3R (1864/437C)	3	-	-	-	-

Note: EIEC - Enteroinvasive E. coli.
 EVEC - Enterovirulent E. coli includes isolates of Enterotoxigenic E. coli, Enteropathogenic E. coli O157 serotypes and non-O157 serotypes.
 Non-EVEC - E. coli isolates not associated with disease.
 The two competitor organisms which are weakly detected by the ompA probes are both Escherichia fergusonii.

25

Equivalents

Those skilled in the art will recognize or be able to ascertain, using no more than routine experimentation, many equivalents to the specific
05 embodiments of the invention described herein. Such equivalents are intended to be encompassed within the scope of this invention.

-88-

SEQUENCE LISTING

(1) GENERAL INFORMATION:

- (i) APPLICANT: Parodos, Kyriaki
McCarty, Janice
- (ii) TITLE OF INVENTION: Nucleic Acid Probes for the Detection of
Shigella
- (iii) NUMBER OF SEQUENCES: 30
- (iv) CORRESPONDENCE ADDRESS:
 - (A) ADDRESSEE: Amoco Corporation
 - (B) STREET: 200 East Randolph Drive, P.O. Box 87703
 - (C) CITY: Chicago
 - (D) STATE: Illinois
 - (E) COUNTRY: U.S.A.
 - (F) ZIP: 60680
- (v) COMPUTER READABLE FORM:
 - (A) MEDIUM TYPE: Floppy disk
 - (B) COMPUTER: IBM PC compatible
 - (C) OPERATING SYSTEM: PC-DOS/MS-DOS
 - (D) SOFTWARE: PatentIn Release #1.0, Version #1.25
- (vi) CURRENT APPLICATION DATA:
 - (A) APPLICATION NUMBER: PCT/US92/06617
 - (B) FILING DATE: 28-JUL-92
- (vii) PRIOR APPLICATION DATA:
 - (A) APPLICATION NUMBER: US 07/738,800
 - (B) FILING DATE: 31-JUL-1991
- (viii) ATTORNEY/AGENT INFORMATION:
 - (A) NAME: Galloway, Norval B.
 - (B) REGISTRATION NUMBER: 33,595
 - (C) REFERENCE/DOCKET NUMBER: GTR90-04 PCT
- (ix) TELECOMMUNICATION INFORMATION:
 - (A) TELEPHONE: 312-856-7180
 - (B) TELEFAX: 312-856-4972

(2) INFORMATION FOR SEQ ID NO:1:

- (i) SEQUENCE CHARACTERISTICS:
 - (A) LENGTH: 164 base pairs
 - (B) TYPE: nucleic acid
 - (C) STRANDEDNESS: double
 - (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:1:

AATCCAACCG CAGTAATAAA CTGAATCCCT CGCATGGCTT GCAGCGCCTC TACTACCGGA	60
TACAGCCTCC ATTCGGTAAC NGCCTCCTTG AGGGCGGATT CCAGCCGTTT ACATTGTGCC	120
TGCCGATCTT CTATTGTACG ACGGTGTTTCG TCAAAGCTA ATTG	164

SUBSTITUTE SHEET

-89-

(2) INFORMATION FOR SEQ ID NO:2:

- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 796 base pairs
 - (B) TYPE: nucleic acid
 - (C) STRANDEDNESS: double
 - (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:2:

CATCAGAAAT CTAAGCAGAA GTCTTTCCGT GTTCTCAGA AATGGGGGCA ACGTGCAAAA	60
CTTGCCCTTG CTGGTGAACA ACGTCTTACA AAGATGGTTC CTGGATGGAT TGACCCTGAG	120
ACTTTTAAAC TCAATGAACA CGCTGAGACT GTGAGATTGA TATTCAAACT GCTGCTAGAT	180
GGTGAAAGTC TGCATAACAT TGCACGTCAC CTTCAAAGCA ACGGTATAAA GTCGTTTAGT	240
CGCCGTAAAG ATGCTAATGG GTTCTCTGTT CACTCTGTAC GCACATTCTA AGGTCAGAGC	300
AACAATAGGC ACGTTACCAG CATCACAACG TAATGACCGC CCCGCTATAC CGAACTACTA	360
CGAAGGTGTT GTAGATATAC CAACGTTCAA TAAAGCTCAA GAGATTCTCG ACAAGAATCG	420
TAAAGGCCGT ACACCTGCAA GTGACAACCC ACTAACGATT AACATCTTCA AAGGTCTGTT	480
TAGGTGTCAG TGTGGGGCTA GTGTCCATCC TACCGGAACA AAGAATAAGT ATGCTGGGGT	540
CTACAGGTGC AATAACCACT TAGACGGTCG CTGTGATGTT CCACCGTTGA AGCGTAAACC	600
GTTTGACCGA TGGATGATTG ATAATTTTCT GGGGATGATT GACGTGGGGA ATGATGGAGA	660
ATCAGAGAGA AAGATTGCAG CGTTACAGCA TGAGGTTGAA ATTGTACAG CCAGAATCAA	720
GAAACGTACC GCCCTACTTC TTGAGATGGA TGATATTGAT GAACTAAAAA TTCAGCTTAA	780
GGAAGTGAAC CAGAAG	796

(2) INFORMATION FOR SEQ ID NO:3:

- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 587 base pairs
 - (B) TYPE: nucleic acid
 - (C) STRANDEDNESS: double
 - (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:3:

GATCTTTCTT CGAAGAAGAG AGCGCACCAA TACCCGCGCC CACGAGAGAG CCCAGACCTG	60
CGCCGATAGC AGATTGCGCT GCTTCGCGTT CGCCGGTGTA AGGGTTAGTT GTGCAGCCAG	120
ATACCGCCAG AGCGCCACTC ACTACGGAGG CAATAAGATA AACACGTTTG TTCATTGTTA	180
ATCCTTCCTA ACCTTTTTAT TCTTGCCAC GGGTCCGTG GCGGGAGATT ATGCCGCGTG	240
AACATGAAGA TGAGGTGTAC TGGCAATAGC GGACACTACC ATTTGTTCTT TTTTAAAGCA	300

SUBSTITUTE SHEET

-90-

GCCATCTGAT GATATTTTTC CCTGAAGGCT GCCGGGGAGA TATTCCCCAG ACGAGAGTGA 360
 CGACGCTGAC GATTCTAGAA AATCTCAATG TATTCCCGTA TTAGTGAGAT GGCTTCATCC 420
 CGGTTATTAA AACGATAGTG GCTCAGGCTC TCATTTTTC ACGTTCCCCA CAAGCTTTCC 480
 ATCGGAGCGT TGTCGTAACA GTTACCTTTA CGCGACATTG ATGTTTTTCAG ACCAGACTGC 540
 TCCTGTATGA CCCGGTAATC GTATGCGCAG TACTGTGAAC CTCGATC 587

(2) INFORMATION FOR SEQ ID NO:4:

- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 786 base pairs
 - (B) TYPE: nucleic acid
 - (C) STRANDEDNESS: double
 - (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:4:

GATCAGAGTG GTGGATTAGC CCGGCAGTGG GCGCTGGCTC CTGAGCGCCA TAAACAGGGC 60
 TTTACCTGTC AGCTCTTTTG TCATGCGCTC TCCCATGGCG TACGACAATT TCGCACGTAT 120
 AAACATCTTT GATGCCAGCG AGGTACAACC ATCCCTCCTG TGTGGCAACA TACGTCAGGT 180
 CCGCCACCCA GACCTGATTT GGTGCTGTAG GAGCGAACGT CTGGTTCAGC AGATTGGCG 240
 CAACTGGCAG ATTGTGGTTC GAGTTCGTAG TCGCTCTGAA CTTGCGTTCT GCTTACAGCG 300
 TAGCTTAGCT CCTTACGAAG ACGTGCCAGT CCGTCACGAC CAACGATGAT GCCATTCTCT 360
 GCCAGCTCCG TCTGGGAGCC GCCGGGTTC CATATGTTTC GCGAGTGGG ATATGTGCCA 420
 CCTTAATCTC CAGTTTTAGC CGCTCATCAC TTTGTTTTCT GTCTGAGGGT TCATGCTGTA 480
 CCCAGTTGTA ATAACCGCTC CTGGATACAC CAAATACTGA CACATCGCTT CAATGGGAAA 540
 TTGTGTGCGC CATTGTTCGA TTAACGCGTA TTTTTCAGCG ACTCCTGTGC AAAATACGCT 600
 GTTGCTTTTT TTAATATATC TCGCTCAAGG CGAGCTTCAT TTAACGCCTT ACGCAGTCGC 660
 AGAATTCAG ATTCAGTTC AGCCACCGTG CGGGAACCGAG GAGTACCGAG CCCTTTTCTG 720
 GCGGCGGTAA CCCATTGTCC TAAAGTGCCT TCAGGAAGAG ATAATCGGGA AGCGCCTTCA 780
 CTGATC 786

(2) INFORMATION FOR SEQ ID NO:5:

- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 1174 base pairs
 - (B) TYPE: nucleic acid
 - (C) STRANDEDNESS: double
 - (D) TOPOLOGY: linear

SUBSTITUTE SHEET

-91-

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:5:

AGTTGGACTA TAGGTGTACT GGCAATACGG ACACACCATT TGTTCCTTTT TTAAGCAGCA	60
TCTGATGATA TTTTCCCTG AAGGCTGCCG GGGAGATATT CCCCAGACGA GAGTGACGAC	120
GCTGACGATT GTAGAAAATC TCAATGTATT CCCGTATTAC TGAGATGGCT TCATCCCGGT	180
TATTAAAACG ATAGTGGCTC AGGCTCTHAT TTTTCAGCGT TCCCCANAAG CTTTCCATCG	240
GAGCGTTGTC GTAACAGTTA CCTTTACGCG ACATTGATGT TTTCAGACCA AACTGCTCCT	300
GTATGACCCG GTAATCGTAT GCGCAGTACT GTGAACCTCG ATCAGAGTGG TGGATTAGCC	360
CGGCAGGNGG GCGCTGGCTC CTGAGCGCCA TAAACAGGGC TTTACCTGTC AGCTCTTTTG	420
TCATGCGCTC TCCCATGCGT AGCCGACAAT TTCGCACGTA TAACATCTTT GATGCCAGCG	480
AGGTACACCA TCCCTCCTGT GTGGCANCAT ACGTCAGGTC CGCCACCCAG ACCTGATTTG	540
GTGCTGTAGG AGTGAACGTC TGGTTCAGCA GATTTGGCGC AACTGGCAGA TTGTGGTTCTG	600
GGTTCGTAGT CGCTCTGAAC TTGCGTTTCT GCTTACAGCG TASCCTAGCT CTTACGAAG	660
ACGTGCCAGT CGGTCACGCC AACGATGATG CCATTCTCTG CCAGCTCCGT CTGGAGCCGC	720
CGGGTTCCAT ATGTTDCGR AGTGCGGATA TGTGCCACCT TAATCTCCAG TTTDAGCCGC	780
TCATCACTTT GTTDTCTGTC TGAGGGTTCA TGCTGTACCC AGTTGTAATA ACCGCTCCTG	840
GATACACCAA ATACCTGACA CATCGCTTSA ATDDDAATT GTTGTCGCCA TTGTTCGATT	900
AACGCGNNNN NNNCAGCGAC TCCTGTGCAA AATACGCTGT TGCTTTTTTTT AATATATCTC	960
GCTCAAGGCG AGCTTCATTT AACGCTTTAC GCAGTTGCAG AATTTAGAT TCCAGTTCAG	1020
CCACCGTGCG GGAACCAGGA GTACCGAGCC CTTTTCTGGC GGCGGTAACC CATTGTCCTA	1080
AAGTGCCTTC AGGAAGAGAT AATCGGGAAG CGCCTTCGCT GATCGAAAGT TGATTTTCAA	1140
GAACCGTTCT GACAGCTTCG GCTTTGAACT CTGT	1174

(2) INFORMATION FOR SEQ ID NO:6:

- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 64 base pairs
 - (B) TYPE: nucleic acid
 - (C) STRANDEDNESS: double
 - (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:6:

TGTGGCATCA ACAATGGTGC GACCACCGAG CGAGATGAGG TGTACTGGCA ATAGCGGACA	60
CAAC	64

SUBSTITUTE SHEET

-92-

(2) INFORMATION FOR SEQ ID NO:7:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 1196 base pairs
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: double
- (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:7:

```

ATTGTAGAAA ATCTCAATGT ATTCCCGTAT TACTGAGATG GCTTCATCCC GGTATTATAA      60
ACGATAGTGG CTCAGGCTCT HATTTTTCAG CGTTCCCCAN AAGCTTTCCA TCGGAGCGTT      120
GTCGTAACAG TTACCTTTAC GCGACATTGA TGTTTTCAGA CCAAACGCT CCTGTATGAC      180
CCGGTAATCG TATGCGCAGT ACTGTGAACC TCGATCAGAG TGGTGGATTA GCGGCGCAGG      240
NGGGCGCTGG CTCCTGAGCG CCATAAACAG GGCTTTACCT GTCAGCTCTT TTGTCATGCG      300
CTCTCCCATG CGTAGCCGAC AATTTGCGAC GTATAACATC TTTGATGCCA GCGAGGTACA      360
CCATCCCTCC TGTGTGGCAN CATACGTCAG GTCGCGCCACC CAGACCTGAT TTGGTGCTGT      420
AGGAGTGAAC GTCTGGTTCA GCAGATTGGG CGGAGCTGGC AGATTGTGGT TCGGGTTTCGT      480
AGTCGCTCTG AACTTGCGTT TCTGCTTACA GCGTASCTTA GCTCCTTACG AAGACGTGCC      540
AGTCGGTCAC GCCAACGATG ATGCCATTCT CTGCCAGCTC CGTCTGGAGC CGCCGGGTTT      600
CATATGTTDN GCRAGTGCGG ATATGTGCCA CCTTAATCTC CAGTTTDAGC CGCTCATCAC      660
TTTGTTDTCT GTCTGAGGGT TCATGCTGTA CCCAGTTGTA ATAACCGCTC CTGGATACAC      720
CAAATACCTG ACACATCGCT TNAATGGGAA ATTGTTGTCG CCATTGTTG ATTAACGCGN      780
ATTTTTCAGC GACTCCTGTG CAAAATACGC TGTGCTTTT TTTNATATAT CTCGCTCAAG      840
GCGAGCTTCA TTTAACGCCT TACGCAGTTG CAGAATTCA GATTCCAGTT CAGCCACCGT      900
GCGGGAACCA GGAGTACCGA GCCCTTTTCT GCGGCGGTA ACCCATGTC CTAAAGTGCC      960
TTCAGGAAGA GNTAATCGGG AAGCGCCTTC GCTGATCGAA AGTTGATTTT CAAGAACCGT     1020
TCTGACAGCT TCGGCTTTGA ACTCTTTAGA GTAACGTTGG TTTTCTGCTC TCATTATTAG     1080
CTCCTTCTGA TGCCATTCTA TTTCAGGAAG GAGTGTCCGT TAAACTCAGG CTACCTCAAG     1140
ATAAAGTTAT TAATTTGCGA GATCACATCT TCAATAGGTT TCGGTCCAT ATTATC         1196

```

(2) INFORMATION FOR SEQ ID NO:8:

(i) SEQUENCE CHARACTERISTICS:

- (A) LENGTH: 64 base pairs
- (B) TYPE: nucleic acid
- (C) STRANDEDNESS: double
- (D) TOPOLOGY: linear

-93-

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:8:

AAAGCACAGA TTTTATAGCT AACTCGATGC TGGTGTGAGG TGTACTGGCA ATAGCGGACA 60
CTAC 64

(2) INFORMATION FOR SEQ ID NO:9:

- (i) SEQUENCE CHARACTERISTICS:
 (A) LENGTH: 1188 base pairs
 (B) TYPE: nucleic acid
 (C) STRANDEDNESS: double
 (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:9:

ATTGTAGAAA ATCTCAATGT ATTCCCGTAT TACTGAGATG GCTTCATCCC GGTTATTAAA 60
 ACGATAGTGG CTCAGGCTCT HATTTTTCAG CGTTCCCCAN AAGCTTTCCA TCGGAGCGTT 120
 GTCGTAACAG TTACCTTTAC GCGACATTRA TGTTTTTCAGA CCAGACTRCT CCTGTATGAC 180
 CCGGTAATCG TATGCGCAGT ACTGTGAACC TCGATCAGAG TGGTGGATTA GCCCGGCAGG 240
 NGGGCGCTGG CTCCTGAGCG CCATAAACAG GGCAAATCCT GTCAGCTCTW TTGTCATGCG 300
 CTCTCCCATG CGTAGCCGAC TAATTTTCGCA CGTATAACAT CTTTGATGCC AGCGAGGTAC 360
 ACCATCCCTC CTGTGTGGCA NCATACGTCA GGTCCGCCAC CCAGACCTGA TTTGGTGCTG 420
 TAGGAGCGAA CGTCTGGTTC AGCAGATTTC GCGCAACTGG CAGATTGTGG TTCGAGTTCG 480
 TAGTCGCTCT GAACTTGCGT TTCTGCTTAC AGTGTAGCCT TAGCTCCTTA CGAAGACGTG 540
 CCAGTCGGTC ACGCCAACGA TGATGCCATT CTCTGCCAGC TCCGTCTGGA GCCGCCGGGT 600
 TCCATATGTT NCGCRAGTGC GGATATGTGC CACCTTAATC TCCAGTTTDA GCCGCTCATC 660
 ACTTTGTTDT CTGTCTGAGG GTTCATGCTG TACCCAGTTG TAATAACCGC TCCTGGATAC 720
 ACCAAATACC TGACACATCG CT TSAATGGG AAATTGTTRT CGCCATTGTT CGATTAACGC 780
 GACTCCTCTG CAAAATACGC TGTGCTTTT TTNATATAT CTGCTCAAG GCGAGCTTCA 840
 TTTAACGCCT TACGCAGTTG CAGAATTTCA GATTCCAGTT CAGCCACCGT GCGGGAACCA 900
 GGAGTACCGA GCCCTTTTCT GGCGGCGGTA ACCCATTGTC CTAAAGTGCC TTCAGGAAGA 960
 GATAATCGGG AAGCGCCTTC ACTGATCGAA AGTTGATTTT CAGGAACCGT TCTGACAGCT 1020
 TCGGCTTTGA ACTCTTKAGA GTAACGTTGG GTTTTCTGCTC TCATTATTAG CTCCTTCTGA 1080
 TGCCATTCTA TTTCAGGAAG GAGTGTCCTG TAAACTCAGG CTACCTCAGT GTGATCGGCG 1140
 ATAAGCCCAG AACTCCGCTC CCAGACCTCC CTGCCAAAAG CAAAACCG 1188

-94-

(2) INFORMATION FOR SEQ ID NO:10:

- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 630 base pairs
 - (B) TYPE: nucleic acid
 - (C) STRANDEDNESS: double
 - (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:10:

```

ATCTCCTTAT GTTATGGAGA TTATTAAAAA GAATAACATT AGCGCTCTCG AACTGCATCG      60
TGCAATTGTT GAGTTGAGTA AAAATATGAA GTCGATTGAT GATAATGCCA GTAAGAAAAA      120
CGACAAGTCA TCATTGTATG TATCATGGAC TCTGAGTTTT ACTGCTCCAA CAAGTAAAGA      180
AGCTCACGAT GTGTTGTCTG GGTATATTAA TTATGTTTCT TCCCTTGTG TAAGGGATTT      240
GATGGAAGAT ATAAGAAATA AACTAGAAGT TAAACTAAT GTTGAAAAAG AAATTCTTGC      300
ACTGGATGAG ATAAAAATTA GAAACCAGCT GAATGCAGAT ATTCGACNCC TCAATTATTC      360
ACTGGAGGTT GCTAATGCGG CTGGAATAAA AAAACCTGTA TACAGCAATG GTCAGATTAT      420
GAAGGATGAC CCAGATTTTC CTGTGGCTCT CGGTTCTGAT GGTATAGCAA CTAAATTGAA      480
CATCAAAAAA TCAATCAAGG ATGTTTCGGA ATTGAGTGGG GAGTTGCGAA ATCGTCAATA      540
TGTGTGAAT CAATTGGTTG TGGCGAAAGN GGGGGANGNN GANNNNANGC MANNNCAGNA      600
NCAANNGTGC CCAACGNNAC CGNCAGAAA      630

```

(2) INFORMATION FOR SEQ ID NO:11:

- (i) SEQUENCE CHARACTERISTICS:
- (A) LENGTH: 388 base pairs
 - (B) TYPE: nucleic acid
 - (C) STRANDEDNESS: double
 - (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:11:

```

CCACCACATT GATTGTCTGC CTGAAATAAC ACGAAAGGCA CTGCGTGAAC GCTATGTGGA      60
ACAGCTGGTG GCTACAGAGA ACAATGTTTC TGAAGTGAAA GCTGTTACCA GAAAAACACG      120
CAATCCTGAC GCTGTCCAGG CAATCGAAGC ATATCGCGGT TCTCCACAAC TGATGGAAGA      180
ACGCCTGAAT GCGCTGACCG AAAACCAGCG CTGGGTATCT GAAGCAAGAG CTGCGCTGGT      240
GGTGGAAGTG CTGAAGCTGG AAAGCGCCGG TAACCCCGGG CGACTGAAAG CCATTAACTT      300
TCTTGTTGAA AAAGCCCCTA AAGGTGAGCT GCCGGAGCGC CTGCAACAGG CCGCAGTTAA      360
CGCCAATGCA AAACGTGGCG CTAATCGT      388

```


-95-

(2) INFORMATION FOR SEQ ID NO:12:

- (i) SEQUENCE CHARACTERISTICS:
(A) LENGTH: 184 base pairs
(B) TYPE: nucleic acid
(C) STRANDEDNESS: double
(D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:12:

```
CTGACGACCT GGACGTGTAC ACTCGTCTGG GTGGTATGGT TTGGCGTGCA GACACCAAAG      60
CTCACAACAA TGTGACAGGT GAATCTGAGA AAAACCACGA TACCGGCGTT TCTCCGGTAT      120
TCGCAGGTGG TGTGAATGG GCCATCACTC CTGAAATCGC TACCCGTCTG GAATACCAGT      180
GGAC                                          184
```

(2) INFORMATION FOR SEQ ID NO:13:

- (i) SEQUENCE CHARACTERISTICS:
(A) LENGTH: 169 base pairs
(B) TYPE: nucleic acid
(C) STRANDEDNESS: double
(D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:13:

```
CTGACGACCT GGACATCTAC ACTCGTCTGG GTGGCATGGT ATGGCGTGCA GACACTAAAT      60
CCAACGTTTA TGGTAAAAAC CACGACACCG GCGTTTCTCC GGTCTTCGCT GGCGGTGTTG      120
AGTACGCGAT CACTCCTGAA ATCGCTACCC GTCTGGAATA CCAGTGGAC                  169
```

(2) INFORMATION FOR SEQ ID NO:14:

- (i) SEQUENCE CHARACTERISTICS:
(A) LENGTH: 35 base pairs
(B) TYPE: nucleic acid
(C) STRANDEDNESS: single
(D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:14:

```
TTGCAGCGCC TCTACTACCG GATACAGCCT CCATT                                  35
```

(2) INFORMATION FOR SEQ ID NO:15:

- (i) SEQUENCE CHARACTERISTICS:
(A) LENGTH: 35 base pairs
(B) TYPE: nucleic acid
(C) STRANDEDNESS: single
(D) TOPOLOGY: linear

-96-

- (xi) SEQUENCE DESCRIPTION: SEQ ID NO:15:
CCTCCTTCAG GCGGATTCC AGCCGTTAC ATTGT 35
- (2) INFORMATION FOR SEQ ID NO:16:
- (i) SEQUENCE CHARACTERISTICS:
(A) LENGTH: 40 base pairs
(B) TYPE: nucleic acid
(C) STRANDEDNESS: single
(D) TOPOLOGY: linear
- (xi) SEQUENCE DESCRIPTION: SEQ ID NO:16:
CCGATCTTCT ATTGTACGAC GGTGTTTCGTC AAAAGCTAAT 40
- (2) INFORMATION FOR SEQ ID NO:17:
- (i) SEQUENCE CHARACTERISTICS:
(A) LENGTH: 41 base pairs
(B) TYPE: nucleic acid
(C) STRANDEDNESS: single
(D) TOPOLOGY: linear
- (xi) SEQUENCE DESCRIPTION: SEQ ID NO:17:
CTGGTGAACA ACGTCTTACA AAGATGGTTC CTGGATGGAT T 41
- (2) INFORMATION FOR SEQ ID NO:18:
- (i) SEQUENCE CHARACTERISTICS:
(A) LENGTH: 41 base pairs
(B) TYPE: nucleic acid
(C) STRANDEDNESS: single
(D) TOPOLOGY: linear
- (xi) SEQUENCE DESCRIPTION: SEQ ID NO:18:
AGTCTTTCCG TGTCTCTCAG AAATGGGGGC AACGTGCAAA A 41
- (2) INFORMATION FOR SEQ ID NO:19:
- (i) SEQUENCE CHARACTERISTICS:
(A) LENGTH: 35 base pairs
(B) TYPE: nucleic acid
(C) STRANDEDNESS: single
(D) TOPOLOGY: linear
- (xi) SEQUENCE DESCRIPTION: SEQ ID NO:19:
CCACCGTTGA AGCGTAAACC GTTTGACCGA TGGAT 35

-97-

(2) INFORMATION FOR SEQ ID NO:20:

- (i) SEQUENCE CHARACTERISTICS:
 - (A) LENGTH: 36 base pairs
 - (B) TYPE: nucleic acid
 - (C) STRANDEDNESS: single
 - (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:20:

GCTGGGGTCT ACAGGTGCAA TAACCACTTA GACGGT

36

(2) INFORMATION FOR SEQ ID NO:21:

- (i) SEQUENCE CHARACTERISTICS:
 - (A) LENGTH: 48 base pairs
 - (B) TYPE: nucleic acid
 - (C) STRANDEDNESS: single
 - (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:21:

CGATGATGCC ATTCTCTGCC AGCTCCGTCT GGGAGCCGCC GGGTTTCC

48

(2) INFORMATION FOR SEQ ID NO:22:

- (i) SEQUENCE CHARACTERISTICS:
 - (A) LENGTH: 17 base pairs
 - (B) TYPE: nucleic acid
 - (C) STRANDEDNESS: single
 - (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:22:

GGAGCAGTCT GGTCTGA

17

(2) INFORMATION FOR SEQ ID NO:23:

- (i) SEQUENCE CHARACTERISTICS:
 - (A) LENGTH: 37 base pairs
 - (B) TYPE: nucleic acid
 - (C) STRANDEDNESS: single
 - (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:23:

CCTGTGGCTC TCGGTTCTGA TGGTATAGCA ACTAAAT

37

-98-

(2) INFORMATION FOR SEQ ID NO:24:

- (i) SEQUENCE CHARACTERISTICS:
(A) LENGTH: 37 base pairs
(B) TYPE: nucleic acid
(C) STRANDEDNESS: single
(D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:24:

CAAGGATGTT TCGGAATTGA GTGGGGAGTT GCGAAAT

37

(2) INFORMATION FOR SEQ ID NO:25:

- (i) SEQUENCE CHARACTERISTICS:
(A) LENGTH: 35 base pairs
(B) TYPE: nucleic acid
(C) STRANDEDNESS: single
(D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:25:

CAGGCAATCG AAGCATATCG CGGTTCTCCA CAACT

35

(2) INFORMATION FOR SEQ ID NO:26:

- (i) SEQUENCE CHARACTERISTICS:
(A) LENGTH: 35 base pairs
(B) TYPE: nucleic acid
(C) STRANDEDNESS: single
(D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:26:

TGAATGCGCT GACCGAAAAC CAGCGCTGGG TATCT

35

(2) INFORMATION FOR SEQ ID NO:27:

- (i) SEQUENCE CHARACTERISTICS:
(A) LENGTH: 35 base pairs
(B) TYPE: nucleic acid
(C) STRANDEDNESS: single
(D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:27:

GTGATGGCCC ATTCAACACC ACCTGCGAAT ACCGG

35

-99-

(2) INFORMATION FOR SEQ ID NO:28:

- (i) SEQUENCE CHARACTERISTICS:
 - (A) LENGTH: 35 base pairs
 - (B) TYPE: nucleic acid
 - (C) STRANDEDNESS: single
 - (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:28:

CTCAGATTCA CCTGTCACAT TGTGTGAGC TTTGG

35

(2) INFORMATION FOR SEQ ID NO:29:

- (i) SEQUENCE CHARACTERISTICS:
 - (A) LENGTH: 41 base pairs
 - (B) TYPE: nucleic acid
 - (C) STRANDEDNESS: single
 - (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:29:

AATCCATCCA GGAACCATCT TTGTAAGACG TTGTTACCA G

41

(2) INFORMATION FOR SEQ ID NO:30:

- (i) SEQUENCE CHARACTERISTICS:
 - (A) LENGTH: 41 base pairs
 - (B) TYPE: nucleic acid
 - (C) STRANDEDNESS: single
 - (D) TOPOLOGY: linear

(xi) SEQUENCE DESCRIPTION: SEQ ID NO:30:

TTTTCACGT TGCCCCATT TCTGAGAAAC ACGGAAAGAC T

41

-100-

CLAIMS

1. An exclusively chromosomal *Shigella* specific
fragment identified by subtractive
hybridization using *S. flexneri* as "target" DNA
05 and a complex DNA competitor mix as competitor
DNA.
2. A *Shigella* specific fragment identified by
subtractive hybridization using *S. sonnei* as
"target" DNA and a mixture of non-pathogenic *E.*
10 coli YMC DNA and pBR322 vector DNA as
competitor DNA.
3. A *Shigella* specific fragment derived from a
chromosomal sequence of *Shigella sonnei*, the
Shigella specific fragment selected from the
15 group consisting of NT-6 (nucleotides 1-124 of
SEQ ID NO:1), NT11-2 (SEQ ID NO:2), NT14 (SEQ
ID NO:4) and NT15 (SEQ ID NO:3).
4. A *Shigella* specific fragment derived from a
chromosomal sequence of *Shigella flexneri*, the
20 Shigella specific fragment selected from the
group consisting of NT18-1a (SEQ ID NO:10) and
NT19-2 (SEQ ID NO:11).
5. A probe for the detection of a member of the
genus *Shigella* wherein the probe is capable of
25 detecting at least one untyped isolate or at
least one isolate of a serotype of each of the

-101-

four species of the genus *Shigella* selected from the group consisting of:

- a) (SEQ ID NO:1);
- b) probes derived from SEQ ID NO:1;
- 05 c) NT11-2 (SEQ ID NO:2);
- d) probes derived from NT11-2; and
- e) homologues of (a) through (d).

6. A probe of Claim 5 for the detection of a member of the genus *Shigella* wherein the probe
- 10 is capable of detecting at least one untyped isolate or at least one isolate of a serotype of each of the four species of the genus *Shigella*, wherein the probe is selected from the group consisting of:

- 15 a) 1500 (SEQ ID NO:14),
5'-TTGCAGCGCCTCTACTACCGGATACAGCCTCCATT-3';
- b) 1501 (SEQ ID NO:15),
5'-CCTCCTTCAGGGCGGATTCCAGCCGTTACATTGT-3';
- c) 1911 (SEQ ID NO:16),
20 5'-CCGATCTTCTATTGTACGACGGTGTTTCGTCAAA-
AGCTAAT-3';
- d) 1682 (SEQ ID NO:17),
5'-CTGGTGAACAACGTCTTACAAAGATGGTTCCTG-
GATGGATT-3';
- 25 e) 1683 (SEQ ID NO:18),
5'-AGTCTTTCCGTGTTTCTCAGAAATGGGGGCAAC-
GTGCAAAA-3';
- f) 1708 (SEQ ID NO:19),
30 5'-CCACCGTTGAAGCGTAAACCGTTTGACCGATGG-
AT-3'; and
- g) complements of (a) through (f).

-102-

7. A probe capable of detecting at least one isolate of each of serotypes 1-5 of *Shigella flexneri* selected from the group consisting of:
- a) NT18-1a (SEQ ID NO:10);
 - 05 b) probes derived from NT18-1a;
 - c) homologues of (a) and (b);
 - d) NT19-2 (SEQ ID NO; 11);
 - e) probes derived from NT19-2; and
 - f) homologues of (d) and (e).
- 10 8. A probe of Claim 7 capable of detecting of at least one isolate of each of serotypes 1-5 of *Shigella flexneri* wherein the probe is selected from the group consisting of:
- a) 1712 (SEQ ID NO:23),
15 5'-CCTGTGGCTCTCGGTTCTGATGGTATAGCAACT-
AAAT-3';
 - b) 1713 (SEQ ID NO: 24),
5'-CAAGGATGTTTCGGAATTGAGTGGGGAGTTGCG-
AAAT-3'; and
 - 20 c) complements of (a) and (b).
9. A probe of Claim 7 (d) through (f) which further detects a member of the genus *Shigella* selected from the group consisting of *S. boydii* serotype 5, *S. boydii* serotype 7, *S. boydii* serotype 9, *S. boydii* serotype 11, *S. boydii* serotype 16, *S. boydii* serotype 17, and *S. sonnei* serotype 1.
- 25

-103-

10. A probe of Claim 9 selected from the group consisting of:
- a) 1684 (SEQ ID NO:25),
5'-CAGGCAATCGAAGCATATCGCGGTTCTCCACAACT-3';
 - 05 b) 1685 (SEQ ID NO: 26),
5'-TGAATGCGCTGACCGAAAACCAGCGCTGGGTATCT-3';
and
 - c) complements of (a) and (b).
11. An oligonucleotide probe for the detection of
10 at least one isolate of *Shigella* selected from the group consisting of *S. dysenteriae* serotypes 1 and 2, and *S. boydii* serotypes 5, 7, 9, 11, 12, 15, 16, and 17, the
oligonucleotide probe comprising a sequence
15 selected from the group consisting of:
- a) 1707 (SEQ ID NO:28),
5'-CTCAGATTACCTGTCACATTGTTGTGAGCTTTGG-3';
 - b) 1706 (SEQ ID NO:27);
5'-GTGATGGCCCATTC AACACCACCTGCGAATACCGG-3';
 - 20 c) complements of (a) and (b); and
 - d) homologues of (a) through (c).
12. A substantially inclusive oligonucleotide probe
for the detection of a member of the genus
Shigella comprising a sequence selected from
25 the group consisting of:
- a) 437 (SEQ ID NO:21),
5'-CGATGATGCCATTCTCTGCCAGCTCCGTCTGG-
GAGCCGCCGGGTTTCC-3';
 - b) 1864 (SEQ ID NO:22),
30 5'-GGAGCAGTCTGGTCTGA-3';

SUBSTITUTE SHEET

-104-

- c) complements of (a) and (b); and
 - d) homologues of (a) and (b).
13. A substantially inclusive capture/detector probe pair for the detection of a member of the genus *Shigella*, comprising a substantially inclusive capture probe and a detector probe wherein the capture probe is selected from the group consisting of:
- a) 1864 (SEQ ID NO:22); and
 - b) complement of 1864;
- the detector probe having an oligonucleotide sequence derived from the same strand of the *S. sonnei* sequence comprising fragments NT14 (SEQ ID NO:4) and NT15 (SEQ ID NO:3) as the capture probe.
14. A substantially inclusive capture/detector probe pair of Claim 13, the probe pair also being exclusive of exclusivity organisms commonly found in stool and substantially exclusive of non-EIEC Enterobacteriaceae.
15. A capture/detector probe pair of Claim 13 wherein the capture probe is 1864 (SEQ ID NO:22) and the detector probe is the complement of 437.
16. A capture/detector probe pair of Claim 13 wherein the capture probe is the complement of 1864 and the detector probe is probe 437 (SEQ ID NO:21).

17. A substantially inclusive probe set for the detection of a member of the genus *Shigella*, which further detects EIEC, comprising at least three capture/detector probe pairs, the first pair of probes selected from the group consisting of oligonucleotides having sequences derived from a selected strand of SEQ ID NO:1 and homologues thereof, the second pair of probes selected from the group consisting of oligonucleotides having sequences derived from a selected strand of fragment NT19-2 (SEQ ID NO:11) and homologues thereof, and the third pair of probes selected from the group consisting of oligonucleotides having sequences derived from a selected strand of the *S. dysenteriae* ompA sequence which spans nucleotides from position 893 through 1076 (SEQ ID NO:12) and homologues thereof.
18. A substantially inclusive probe set of Claim 17, comprising at least three capture/detector probe pairs, the first pair of probes having oligonucleotide sequences derived from a selected strand of SEQ ID NO:1, the second pair of probes having oligonucleotide sequences derived from a selected strand of fragment NT19-2 (SEQ ID NO: 11), and the third pair of probes having oligonucleotide sequences derived

-106-

from a selected strand of the *S. dysenteriae* ompA sequence which spans nucleotides from position 893 through 1076 (SEQ ID NO:12).

19. A substantially inclusive probe set of Claim
05 18, the probe set also being exclusive of
exclusivity organisms commonly found in stool
and substantially exclusive of non-EIEC
Enterobacteriaceae.
20. A substantially inclusive probe set of Claim 19
10 wherein the first pair of probes is selected
from the group consisting of probes 1911 (SEQ
ID NO:16), 1500 (SEQ ID NO:14) and 1501 (SEQ ID
NO:15), the second pair of probes consists of
15 probes 1684 (SEQ ID NO:25) and 1685 (SEQ ID
NO:26), and the third pair of probes consists
of probes 1706 (SEQ ID NO:27) and 1707 (SEQ ID
NO:28).
21. A substantially inclusive probe set for the
20 detection of a member of the genus *Shigella*,
which further detects EIEC, comprising at least
three capture/detector probe pairs, the first
pair selected from the group consisting of:
a) 1911 (SEQ ID NO:16);
b) 1500 (SEQ ID NO:14); and
25 c) 1501 (SEQ ID NO:15);
the second pair selected from the group
consisting of:
a) 1684 (SEQ ID NO:25); and

-107-

- b) 1685 (SEQ ID NO:26);
the third pair selected from the group
consisting of:
a) 1706 (SEQ ID NO:27); and
05 b) 1707 (SEQ ID NO:28).
22. A substantially inclusive probe set for the
detection of a member of the genus *Shigella*,
which further detects EIEC, comprising at least
two capture/detector probe pairs, the first
10 pair of probes selected from the group
consisting of oligonucleotides having sequences
derived from a selected strand of the sequence
comprising *S. sonnei* fragments NT14 (SEQ ID
NO:4) and NT15 (SEQ ID NO:3) and homologues
15 thereof, and the second pair of probes selected
from the group consisting of oligonucleotides
having sequences derived from a selected strand
of the *S. dysenteriae* ompA sequence which spans
nucleotides from position 893 through 1076 (SEQ
20 ID NO:12) and homologues thereof.
23. A substantially inclusive probe set of Claim
22, comprising at least two capture/detector
probe pairs, the first pair of probes having
oligonucleotide sequences derived from a
25 selected strand of the sequence comprising *S.*
sonnei fragments NT14 (SEQ ID NO:4) and NT15
(SEQ ID NO:3), and the second pair of probes
having oligonucleotide sequences derived from a
selected strand of the *S. dysenteriae* ompA
30 sequence which spans nucleotides from position
893 through 1076 (SEQ ID NO:12).

SUBSTITUTE SHEET

-108-

24. A substantially inclusive probe set of Claim
23, the probe set also being exclusive of
exclusivity organisms commonly found in stool
and substantially exclusive of non-EIEC
05 Enterobacteriaceae.
25. A substantially inclusive probe set of Claim 24
wherein the first pair of probes consists of
probes 1864 (SEQ ID NO:22) and the complement
of 437, and the second pair of probes consists
10 of probes 1706 (SEQ ID NO:27) and 1707 (SEQ ID
NO:28).
26. A substantially inclusive probe set for the
detection of a member of the genus Shigella,
which further detects EIEC, comprising at least
15 two capture/detector probe pairs, wherein the
first pair is 1864 (SEQ ID NO:22) and the
complement of 437 and the second pair is 1706
(SEQ ID NO:27) and 1707 (SEQ ID NO:28).
27. An oligonucleotide probe derived from Shigella
20 specific fragment NT-6 (nucleotides 1-124 of
SEQ ID NO:1) which substantially retains the
inclusivity behavior of NT-6.
28. An oligonucleotide probe of Claim 27 which has
improved exclusivity behavior for non-EIEC
25 Enterobacteriaceae and is exclusive of the
exclusivity organisms commonly found in stool.

29. An oligonucleotide probe of Claim 28 comprising a sequence selected from the group consisting of:
- 05 a) 1911 (SEQ ID NO:16),
5'-CCGATCTTCTATTGTACGACGGTGTTCGTCAAA-
AGCTAAT-3';
- b) 1500 (SEQ ID NO:14),
5'-TTGCAGCGCCTCTACTACCGGATACAGCCTCCATT-3';
- 10 c) 1501 (SEQ ID NO:15),
5'-CCTCCTTCAGGGCGGATTCCAGCCGTTACATTGT-3';
- d) complements of (a) through (c); and
- e) homologues of (a) through (d).
30. An oligonucleotide probe derived from *Shigella* specific fragment NT11-2 (SEQ ID NO:2) which
- 15 substantially retains the inclusivity behavior of NT11-2.
31. An oligonucleotide probe of Claim 30 which substantially retains the exclusivity behavior of NT11-2 (SEQ ID NO:2) towards non-EIEC
- 20 Enterobacteriaceae.

-110-

32. An oligonucleotide probe of Claim 31 comprising a sequence selected from the group consisting of:
- a) 1682 (SEQ ID NO:17),
- 05 5'-CTGGTGAACAACGTCTTACAAAGATGGTTCCTG-
GATGGATT-3';
- b) complement of 1682 (SEQ ID NO:29),
5'-AATCCATCCAGGAACCATCTTTGTAAGACGTTG-
TTCACCAG-3'; and
- 10 c) homologues of (a) and (b).
33. An oligonucleotide probe derived from *Shigella* specific fragment NT11-2 (SEQ ID NO:2) which moderately retains the inclusivity behavior of NT11-2.
- 15 34. An oligonucleotide probe of Claim 33 which substantially retains the exclusivity behavior of NT11-2 (SEQ ID NO:2) towards non-EIEC Enterobacteriaceae.
- 20 35. An oligonucleotide probe of Claim 34 comprising a sequence selected from the group consisting of:
- a) 1683 (SEQ ID NO:18),
5'-AGTCTTTCCGTGTTTCTCAGAAATGGGGGCAAC-
GTGCAAAA-3';
- 25 b) complement of 1683 (SEQ ID NO:30),
5'-TTTTGCACGTTGCCCCATTCTGAGAAACACGG-
AAAGACT-3'; and
- c) homologues of (a) and (b).

-111-

36. An oligonucleotide probe derived from *Shigella* specific fragment NT11-2 (SEQ ID NO:2) which partially retains the inclusivity behavior of NT11-2.
- 05 37. An oligonucleotide probe of Claim 36 which substantially retains the exclusivity behavior of NT11-2 (SEQ ID NO:2) towards non-EIEC Enterobacteriaceae.
- 10 38. An oligonucleotide probe of Claim 37 comprising a sequence selected from the group consisting of:
- a) 1708 (SEQ ID NO:19),
5'-CCACCGTTGAAGCGTAAACCGTTTGACCGATGG-
AT-3'; and
- 15 b) 1709 (SEQ ID NO:20),
5'-GCTGGGGTCTACAGGTGCAATAACCACTTAGAC-
GGT'-3';
- c) complements of (a) and (b); and
- d) homologues of (a) through (c).
- 20 39. An oligonucleotide probe derived from *Shigella* specific fragment NT18-1a (SEQ ID NO:10) which substantially retains the inclusivity behavior of NT18-1a.
- 25 40. An oligonucleotide probe of Claim 39 which substantially retains the exclusivity behavior of NT18-1a (SEQ ID NO:10) towards non-EIEC Enterobacteriaceae.

SUBSTITUTE SHEET

-112-

41. An oligonucleotide probe of Claim 40 comprising a sequence selected from the group consisting of:
- 05 a) 1712 (SEQ ID NO:23),
5'-CCTGTGGCTCTCGGTTCTGATGGTATAGCAACT-
AAAT-3';
- b) 1713 (SEQ ID NO:24),
5'-CAAGGATGTTTCGGAATTGAGTGGGGAGTTGCG-
AAAT-3';
- 10 c) complements of (a) and (b); and
d) homologues of (a) through (c).
42. An oligonucleotide probe derived from *Shigella* specific fragment NT19-2 (SEQ ID NO:11) which substantially retains the inclusivity behavior
- 15 of NT19-2.
43. An oligonucleotide probe of Claim 42 which substantially retains the exclusivity behavior of NT19-2 (SEQ ID NO:11) towards non-EIEC *Enterobacteriaceae*.
- 20 44. An oligonucleotide probe of Claim 43 comprising a sequence selected from the group consisting of:
- a) 1684 (SEQ ID NO:25),
5'-CAGGCAATCGAAGCATATCGCGGTTCTCCACAACT-3';
- 25 b) 1685 (SEQ ID NO:26),
5'-TGAATGCGCTGACCGAAAACCAGCGCTGGGTATCT-3';
- c) complements of (a) and (b); and
d) homologues of (a) through (c).

1/15




FIG. 1

2/15

FIGURE 2

NT6	5'	aaTCCAACCGCAGTAATAAACTGAATCCCTCGCATGGCTTGACAGCGCCTCTACTA	
Probe 1500		TTGCAGCGCCTCTACTA	
NT6 cont'd		CCGGATACAGCCTCCATTCCGGTAACN9CCTTCAGGGCGGATTCCAGCCGTTT	
1500 cont'd		CCGGATACAGCCTCCATT	
Probe 1501		CCTCCTTCAGGGCGGATTCCAGCCGTTT	
NT6 cont'd		ACATTGTGCCCTGCCGATCTTCTATTTGTACGACGGTGTTCGTCAAAAGCTAATTG	3'
1501 cont'd		ACATTGT	
Probe 1911		CCGATCTTCTATTTGTACGACGGTGTTCGTCAAAAGCTAAT	

Sau3A

| | |

SUBSTITUTE SHEET

FIGURE 3A

NT11-2 5' CATCAGAAATCTAAGCAGAAAGTCTTTCCGTGTTTCTCAGAAATGGGGCAACGTGC
Probe 1683 AGTCTTTCCGTGTTTCTCAGAAATGGGGCAACGTGC

NT11-2 cont'd AAAACTTGCCCTTGCTGGTGAAACAACGTCTTACAAAGATGGTTCCTGGATGGATTG
1683 cont'd AAAA
Probe 1682 CTGGTGAAACAACGTCTTACAAAGATGGTTCCTGGATGGATT

NT11-2 cont'd ACCCTGAGACTTTTAAACTCAATGAACACGCTGA/ACTGTGAGATTGATATTCAAA

NT11-2 cont'd CTGCTGCTAGATGGTGAAAGTCTGCATAACATTGCACGTCACCTTCAAAGCAACGG

NT11-2 cont'd TATAAAGTCGTTTAGTCGCCCGTAAAGATGCTAATGGGTtCTCTGTCTCACTCTGTAC

NT11-2 cont'd GCACATTCTAAGGTCAGAGCAACAATAGGCACGTTACCAGCATCACAAACGTAATGA

NT11-2 cont'd CCGCCCCGCTATACCGAACTACTACGAAGGTGTTGTAGATATACCAACGTTCAATA

NT11-2 cont'd AAGCTCAAGAGATTCTCGACAAGAATCGTAAAGGCCGTACACCTGCAAGTGACAAC

SUBSTITUTE SHEET

4/15

FIGURE 3B

NT11-2 cont'd	CCACTAACGATTAAACATCTTCAAAGGtCTGTTAGGTGTCA GTGTGGGGCTAGTGT
NT11-2 cont'd	CCATCCTACCGGAACAAAGAAATAAGTATGCTGGGGTCTACAGGTGCAATAACCACT
Probe 1709	GCTGGGGTCTACAGGTGCAATAACCACT
NT11-2 cont'd	TAGACGTCGCTGTGATGTTCCACCGTTGAAGCGTAAACCGTTTGACCGATGGATG
1709 cont'd	TAGACGGT
Probe 1708	CCACCGTTGAAGCGTAAACCGTTTGACCGATGGAT
NT11-2 cont'd	ATTGATAATTTTCTGGGGATGATTGACGTGGGGGAATGATGGAGAATCAGAGAGAA
NT11-2 cont'd	AGATTGCAGCGTTACAGCATGAGGTTGAAATTGTCA CAGCCAGAAATCAAGAAACG
NT11-2 cont'd	TACCGCCCTACTTCTTGAGATGGATGATATTGATGA ACTAAAAATTCAGCTTAAG
NT11-2 cont'd	3' GAACTGAAC CAGAAG

SUBSTITUTE SHEET

5/15

FIGURE 4A

5'
 NT15 GATCTTCTTCGAAGAAGAGAGCGCACCAATACCCGCGCCACGAGAGAGCCCAG
 NT15 ACCTGCGCCGATAGCAGATTTCCTGCTTCGCGTTTCGCGGTGTAAGGGTTAGTT
 NT15 GTGCAGCCAGATACCGCCAGAGCGCCACTCACTACGGAGGCAATAAGATAAACAC
 NT15 GTTTGTTCAATTGTTAATCCTTCCTAACCTTTTATTCTTTGCCACGGGTTCCGTG
 E.c. 2 TGTGG
 S.f. AAAGC
 NT15 GCGGGAGATTATGCCGCGTGAACATGAAGATGAGGTGTACTGGCAATAGCGGACA
 E.c. 1 AGTGGACTATAGTgTACTgGCAATA-CgGACA
 E.c. 2 CATCAACAATGGTGCAGACCCGAGCGAGATGAGGTGTACTGGCAATAGCGGACA
 S.f. ACAGATTTTATAGCTAACTCGATGCTGGTGTGAGGTGTACTGGCAATAGCGGACA

 Left end of repeat.
 NT15 CTACCATTTGTTCTTTTTTAAGCAGCCATCTGATGATATTTTCCCTGAAGGCT
 E.c. 1 C-ACCATTTGTTctttttTTAAGCAG-CATCTGATGATAtTTTCCCTGAAGGCT
 E.c. 2 CAAC-----
 S.f. CTAC-----

SUBSTITUTE SHEET

6/15

FIGURE 4B

NT15
E.C. 1
E.C. 2
S.f.
GCCGGGAGATATCCCCAGACGAGTGACGACGCTGACGATTCTAGAAAAATCTC
GCCgGGAGATATTCCCCAGACGAGTGACGACGCTGACGATTGTAGAAAAATCTC
-----ATTGTAGAAAAATCTC
-----ATTGTAGAAAAATCTC

NT15
E.C. 1
E.C. 2
S.f.
AATGTATTCCCGTATTACTGAGATGGCTTTCATCCCGGTTATTAAACGATAGTGGC
AATGTATTCCCGTATTACTGAGATGGCTTTCATCCCGGTTATTAAACGATAGTGGC
AATGTATTCCCGTATTACTGAGATGGCTTTCATCCCGGTTATTAAACGATAGTGGC
AATGTATTCCCGTATTACTGAGATGGCTTTCATCCCGGTTATTAAACGATAGTGGC

NT15
E.C. 1
E.C. 2
S.f.
TCAGGCTCTCATTTTTCAGCGTTCCCCACAAAGCTTTCCATCGGAGCGTTGTCGTAA
TCAGGCTCTcHATTTTTCAGCGTTCCCCCANAAAGCTTTCCATCGGAGCGTTGTCGTAA
TCAGGCTCTcHATTTTTCAGCGTTCCCCCANAAAGCTTTCCATCGGAGCGTTGTCGTAA
TCAGGCTCTcHATTTTTCAGCGTTCCCCCANAAAGCTTTCCATCGGAGCGTTGTCGTAA

NT15
E.C. 1
E.C. 2
S.f.
CAGTTACCTTTACGGACATTGATGTTTTCAGACCCAGACTGCTCCTGTATGACCCGG
CAGTTACCTTTACGGACATTGATGTTTTCAGACCCAAACTGCTCCTGTATGACCCGG
CAGTTACCTTTACGGACATTGATGTTTTCAGACCCAAACTGCTCCTGTATGACCCGG
CAGTTACCTTTACGGACATTGATGTTTTCAGACCCAGACTRCTCCTGTATGACCCGG
3' AGTCTGGTCTGACGAGG 5'

Probe 1864



FIGURE 4C

NT15	TAAATCGTATGCGCAGTACTGTGAACCTCGATC	
NT14		GATCAGAGTGGTGGATTAGCCCGGCAG-
E.C. 1	TAAATCGTATGCGCAGTACTGTGAACCTCGATCAGAGTGGTGGATTAGCCCGGCAGG	
E.C. 2	TAAATCGTATGCGCAGTACTGTGAACCTCGATCAGAGTGGTGGATTAGCCCGGCAGG	
S.f.	TAAATCGTATGCGCAGTACTGTGAACCTCGATCAGAGTGGTGGATTAGCCCGGCAGG	
NT14	TGGGCGCTGGCTCCTGAGCGCCATAAACAGGGCTTTACCTGTCAAGCTCTTTTGTCA	
E.C. 1	NGGGCGCTGGCTCCTGAGCGCCATAAACAGGGCTTTACCTGTCAAGCTCTTTTGTCA	
E.C. 2	NGGGCGCTGGCTCCTGAGCGCCATAAACAGGGCTTTACCTGTCAAGCTCTTTTGTCA	
S.f.	NGGGCGCTGGCTCCTGAGCGCCATAAACAGGGCAATCCTGTCAAGCTCTWTTGTCA	
NT14	TGCGCTCTCCCATGGCGTA--CGAC-AATTTGCGACGTATAAACATCTTTGATGCC	
E.C. 1	TGCGCTCTCCCATG-CGTAGCCGAC-AATTTGCGACGTATAA-CATCTTTGATGCC	
E.C. 2	TGCGCTCTCCCATG-CGTAGCCGAC-AATTTGCGACGTATAA-CATCTTTGATGCC	
S.f.	TGCGCTCTCCCATG-CGTAGCCGACTAATTTGCGACGTATAA-CATCTTTGATGCC	
NT14	AGCGAGGTACAACCATCCCTCTGTGTGGCAACATACGTCAAGTCCGCCACCCAGA	
E.C. 1	AGCGAGGTACA-CCATCCCTCCTGTGTGGCANCATACGTCAAGTCCGCCACCCAGA	
E.C. 2	AGCGAGGTACA-CCATCCCTCCTGTGTGGCANCATACGTCAAGTCCGCCACCCAGA	
S.f.	AGCGAGGTACA-CCATCCCTCCTGTGTGGCANCATACGTCAAGTCCGCCACCCAGA	

8/15

FIGURE 4D

NT14		CCTGATTGGTGCTGTAGGAGCGAACGTCTGTGTTTCAGCAGATTGGCGCAACTGGC
E.c. 1		CCTGaTTGGTGCTGTAGGAGTGAAACGTCTGTGTTTCAGCAGATTGGCGCAACTGGC
E.c. 2		CCTGaTTGGTGCTGTAGGAGTGAAACGTCTGTGTTTCAGCAGATTGGCGCAACTGGC
S.f.		CCTGaTTGGTGCTGTAGGAGCGAACGTCTGTGTTTCAGCAGATTGGCGCAACTGGC
NT14		AGATTGTGGTTCGAGTTCGTAGTCGCTCTGAACCTTGCGTT-CTGCTTACAGCGTAG
E.c. 1		AGATTGTGGTTCGGGTTTCGTAGTCGCTCTGAaCTTGCGTTTCTGCTTACAGCGTAS
E.c. 2		AGATTGTGGTTCGGGTTTCGTAGTCGCTCTGAaCTTGCGTTTCTGCTTACAGCGTAS
S.f.		AGATTGTGGTTCGAGTTCGTAGTCGCTCTGAaCTTGCGTTTCTGCTTACAGTGTAG
NT14		-CTTAGCTCCTTACGAAGACGTGCCAGTCGGTCACGACCAACGATGCCATTCT
E.c. 1		-CTTAGCTCCTTACGAAGACGTGCCAGTCGGTCACG-CCAACGATGATGCCATTCT
E.c. 2		-CTTAGCTCCTTACGAAGACGTGCCAGTCGGTCACG-CCAACGATGATGCCATTCT
S.f.		CCTTAGCTCCTTACGAAGACGTGCCAGTCGGTCACG-CCAACGATGATGCCATTCT
Probe 437	5'	CGATGATGCCATTCT
NT14		CTGCCAGCTCCGTCTGG-AGCCGCCGGGTT-CCATATGTTTCGCGAGTCCGGATAT
E.c. 1		CTGCCAGCTCCGTCTGG-AGCCGCCGGGTT-CCATATGTTTCGCGAGTCCGGATAT
E.c. 2		CTGCCAGCTCCGTCTGG-AGCCGCCGGGTT-CCATATGTTDNGCRAGTCCGGATAT
S.f.		CTGCCAGCTCCGTCTGG-AGCCGCCGGGTT-CCATATGTTNCGCRAGTCCGGATAT
437 cont'd	3'	CTGCCAGCTCCGTCTGGGAGCCGCCGGGTTTCC

9/15

FIGURE 4E

NT14		GTGCCACCTTAATCTCCAGTTTTAGCCGCTCATCACTTTGTTTCTGTCTGAGGGT
E.C.	1	GTGCCACCTTAATCTCCAGTTTDAAGCCGCTCATCACTTTGTTTCTGTCTGAGGGT
E.C.	2	GTGCCACCTTAATCTCCAGTTTDAAGCCGCTCATCACTTTGTTTCTGTCTGAGGGT
S.f.		GTGCCACCTTAATCTCCAGTTTDAAGCCGCTCATCACTTTGTTTCTGTCTGAGGGT
NT14		TCATGCTGTACCCAGTTGTAATAACCGCTCCTGGATACACCAAATAC-TGACACAT
E.C.	1	tCATGCTGTACCCAGTTGTAATAACCGCTCCTGGATACACCAAATACCTGACACAT
E.C.	2	TCATGCTGTACCCAGTTGTAATAACCGCTCCTGGATACACCAAATACCTGACACAT
S.f.		TCATGCTGTACCCAGTTGTAATAACCGCTCCTGGATACACCAAATACCTGACACAT
NT14		CGCTTCAATGGGAAATTGTTGTCGCCATTGTTTCGATTAAACGCGTATTTTTCAGCGA
E.C.	1	CGCTTSAATDDDAATTGTTGTCGCCATTGTTTCGATTAAACGCG-----CAGCGA
E.C.	2	CGCTTNAATGGGAAATTGTTGTCGCCATTGTTTCGATTAAACGCGnATTTTTCAGCGA
S.f.		CGCTTSAATGGGAAATTGTTTCGCCATTGTTTCGATTAAACGCG-----A
NT14		CTCCTGTGC AAAAATACGCTGTTGCTTTTAAATATATCTCGCTCAAGCGGAGCTT
E.C.	1	CTCCTGTGC AAAAATACGCTGTTGCTTTTAAATATATCTCGCTCaAGCGGAGCTT
E.C.	2	CTCCTGTGC AAAAATACGCTGTTGCTTTTAAATATATCTCGCTCAAGCGGAGCTT
S.f.		CTCCTGTGC AAAAATACGCTGTTGCTTTTAAATATATCTCGCTCaAGCGGAGCTT
NT14		CATTTAACGCCCTTACGCAGTGCAGAAATTTCAGATTCCAGTTCAGCCACCGTGCGG
E.C.	1	CATTTAACGCCCTTACGCAGTGCAGAAATTTCAGATTCCAGTTCAGCCACCGTGCGG
E.C.	2	CATTTAACGCCCTTACGCAGTGCAGAAATTTCAGATTCCAGTTCAGCCACCGTGCGG
S.f.		CATTTAACGCCCTTACGCAGTGCAGAAATTTCAGATTCCAGTTCAGCCACCGTGCGG

10/15

FIGURE 4F

NT14		GAACCAGGAGTACCGAGCCCTTTTCTGGCGGGTAACCCATTGTCTCCTAAAGTGCC
E.C.	1	GAaCCAGGAGTACCGAGCCCTTTTCTGGCGGGTAaCCCATTTGTCTCCTAAAGTGCC
E.C.	2	GAaCCAGGAGTACCGAGCCCTTTTCTGGCGGGTAaCCCATTTGTCTCCTAAAGTGCC
S.f.		GAaCCAGGAGTACCGAGCCCTTTTCTGGCGGGTAaCCCATTTGTCTCCTAAAGTGCC
NT14		TTCAGGAAGAGATAATCGGGAAGCGCCTTCACTGATC
E.C.	1	TTCAGGAAGAGaTAATCGGGAAGCGCCTTTCGCTGATCGAAAGTTGATTTTCAAGAA
E.C.	2	TTCAGGAAGAGnTAATCGGGAAGCGCCTTTCGCTGATCGAAAGTTGATTTTCAAGAA
S.f.		TTCAGGAAGAGaTAATCGGGAAGCGCCTTTCGCTGATCGAAAGTTGATTTTCAAGAA
E.C.	1	CCGTTCTGACAGCTTCGGCCTTTGAACTCTgT
E.C.	2	CCGTTCTGACAGCTTCGGCCTTTGAACTCTTAGAGTAACGTTGGTTTCTGCTC
S.f.		CCGTTCTGaCAGCTTCGGCCTTTGAACTCTTKAGAGtAACGTTGGGTTTCTGCTC
E.C.	2	ATTATTAGCTCCTTCTGATGCCATTCTATTTCAGGAAGGAGTGTCCTTAAACTCA
S.f.		ATTATTAGCTCCTTCTGATGCCATTCTATTTCAGGAAGGAGTGTCCTTAAACTCA
E.C.	2	GGCTACCTCAAGATAAAGTTATTAAATTTTCGAAGATCACATCTTCAATAGGTTTGGC
S.f.		GGCTACCTCAGTGTGATCGGCGATAAGCCCAgAACTCCGCTCCCAGACCTCCCTGC

		Right end of repeat.
E.C.	2	GTCCATATTATC
S.f.		CAAAAGCAAAACCG
		3'

11/15

FIGURE 5A

NT18-1a	5'	ATCTCCTTATGTTATGGAGATTATTAAAAAGAAATAACATTAGCGCTCTCGAACT
NT18-1a	cont'd	GCATCGTGCAATTGTTGAGTTGAGTAAAAAATATGAAGTCGATTGATGATAATGC
NT18-1a	cont'd	CAGTAAGAAAAACGACAAGTCATCATTTGTATGTATCATGGACTCTGAGTTTAC
NT18-1a	cont'd	TGCTCCAACAAGTAAAGAAGCTCACGATGTGTTGTCTGGGTATATTAAATTATGT
NT18-1a	cont'd	TTCTTCCCTTGTGTAAGGGATTGTGATGGAAAGATATAAGAAATAAACTAGAAAGT
NT18-1a	cont'd	TAAAACTAATGTTGAAAAAAGAAATTCTTGCACTGGATGAGATAAAAAATTAGAAA
NT18-1a	cont'd	CCAGCTGAATGCAGATATTTCGACNCCTCAATTATTCACTGGAGGTTGCTAATGC
NT18-1a	cont'd	GGCTGGAATAAAAAACCTGTATACAGCAATGGTCAGATTATGAAGGATGACCC

SUBSTITUTE SHEET

12/15

FIGURE 5B

NT18-1a cont'd Probe 1712	AGATTTTCCTGTGGCTCTCGGTTCTGATGGTATAGCAACTAAATTGAACATCAA CCTGTGGCTCTCGGTTCTGATGGTATAGCAACTAAAT
NT18-1a cont'd Probe 1713	AAAATCAATCAAGGATGTTTCGGAATTGAGTGGGAGTTGCCGAAATCGTCAATA CAAGGATGTTTCGGAATTGAGTGGGAGTTGCCGAAAT
NT18-1a cont'd	TGTTGTGAATCAATTGGTTGTGGCGAAAGNCGGGGANGNNGANNNNANGCMANN
NT18-1a cont'd	NCAGNANCAANNGTGCCCAACGNNACCGGNCagaaa 3'

13/15

FIGURE 6A

NT19-2	5'	CCACCACATTGATTGCTGCTGAAATAACACGAAAGGCACtGCCGTGAACGCTA
NT19-2 cont'd		TGTGGAACAGCTGGTGGCTACAGAGAAACAATGTTTCTGAAGTGAAAGCTGTTAC
NT19-2 cont'd		CAGAAAAACACGCAATCCTGACGCTGTCCAGGCAATCGAAGCATATCGCGGTTC
Probe 1684		CAGGCAATCGAAGCATATCGCGGTTC
NT19-2 cont'd		TCCACAACTGATGGAAGAACGCCCTGAATGCCGCTGACCGAAAAACAGCGCTGGGT
1684 cont'd		TCCACAACT
Probe 1685		TGAATGCCGCTGACCGAAAAACAGCGCTGGGT
NT19-2 cont'd		ATCTGAAGCAAGAGCTGCCGCTGGTGGTGAAGTGCTGAAGCTGGAAAGCGCCGG
1685 cont'd		ATCT

14/15

FIGURE 6B

NT19-2 cont'd	TAACCCCGGCGACTGAAAGCCATTAACTTCTTGTGTGAAAAAGCCCTAAAGG
NT19-2 cont'd	TGAGCTGCCGGAGCGCCTGCAACAGGCCGAGTTAACGCCAATGC AAAACGTGG
NT19-2 cont'd	CGCTAATCGT 3'

FIGURE 7

S.d. ompA	5'	CTGACGACCTGGACGTGTACACTCGTCTGGTGGTATGGTTTGGCGTGCAGA
E.c. ompA		CTGACGACCTGGACATCTACACTCGTCTGGTGGCATGGTATGGCGTGCAGA
S.d. ompA cont'd		CACCAAAGCTCACAAATGTGACAGGTGAATCTGAGAAAAACCGATACC
E.c. ompA cont'd		CACTAAATCC-----AACGTTTATGGT-----AAAAACCGACACC
Probe 1707	3'	GGTTTCGAGTGTGTTACACTGTCCACTTAGACTC 5'
S.d. ompA cont'd		GGCGTTTCTCCGGTATTCCGACAGGTGGTGTGAATGGGCCATCACTCCTGAAA
E.c. ompA cont'd		GGCGTTTCTCCGGTCTTCGCTGGCGGTGTGAGTACGGCATCACTCCTGAAA
Probe 1706		3' GGCCATAAGCGTCCACCACTTACCCGGTAGTG 5'
S.d. ompA cont'd		TCGCTACCCGCTCTGGAATACCAAGTGGAC
E.c. ompA cont'd		TCGCTACCCGCTCTGGAATACCAAGTGGAC

3'

SUBSTITUTE SHEET

PCTWORLD INTELLECTUAL PROPERTY ORGANIZATION
International Bureau

INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

(51) International Patent Classification ⁵ : C12Q 1/68	A3	(11) International Publication Number: WO 93/03187 (43) International Publication Date: 18 February 1993 (18.02.93)
(21) International Application Number: PCT/US92/06617 (22) International Filing Date: 28 July 1992 (28.07.92) (30) Priority data: 738,800 31 July 1991 (31.07.91) US (71) Applicant: AMOCO CORPORATION [US/US]; Mail Code 1907, Patents and Licensing Department, P.O. Box 87703, Chicago, IL 60680-0703 (US). (72) Inventors: PARODOS, Kyriaki ; 30 Royal Crest Drive, Apartment 9, Marlborough, MA 01752 (US). McCARTY, Janice, Marie ; 10 Donna Terrace, Hyde Park, MA 02136 (US).		(74) Agents: LADD, Robert, G. et al.; Amoco Corporation, Patents and Licensing Department, Mail Code 1907, P.O. Box 87703, Chicago, IL 60680-0703 (US). (81) Designated States: JP, European patent (AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LU, MC, NL, SE). Published <i>With international search report.</i> <i>Before the expiration of the time limit for amending the claims and to be republished in the event of the receipt of amendments.</i> (88) Date of publication of the international search report: 1 April 1993 (01.04.93)
(54) Title: NUCLEIC ACID PROBES FOR THE DETECTION OF SHIGELLA (57) Abstract The invention relates to methods of detection of bacteria of the genus <i>Shigella</i> and/or Enteroinvasive <i>E. coli</i> (EIEC) by use of a set of nucleic acid probes. The invention further relates to a set of <i>Shigella</i> specific chromosomal sequences and fragments and to probes derived from the <i>Shigella</i> specific fragments. Additionally, probes were derived from a sequence from the <i>Shigella</i> <i>ompA</i> gene. In particular, a series of probes, each approximately 40 nucleotides in length, were designed having specificity for <i>Shigella</i> or for <i>Shigella</i> and Enteroinvasive <i>E. coli</i> , and having utility in nonisotopic test formats which require amplification to achieve high sensitivity. Specific hybridization probe sets which are capable of detecting substantially all clinically significant serotypes of <i>Shigella</i> , as well as enteroinvasive strains of <i>E. coli</i> , are disclosed.		

FOR THE PURPOSES OF INFORMATION ONLY

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AT	Austria	FI	Finland	MN	Mongolia
AU	Australia	FR	France	MR	Mauritania
BB	Barbados	GA	Gabon	MW	Malawi
BE	Belgium	GB	United Kingdom	NL	Netherlands
BF	Burkina Faso	GN	Guinea	NO	Norway
BG	Bulgaria	GR	Greece	NZ	New Zealand
BJ	Benin	HU	Hungary	PL	Poland
BR	Brazil	IE	Ireland	PT	Portugal
CA	Canada	IT	Italy	RO	Romania
CF	Central African Republic	JP	Japan	RU	Russian Federation
CG	Congo	KP	Democratic People's Republic of Korea	SD	Sudan
CH	Switzerland	KR	Republic of Korea	SE	Sweden
CI	Côte d'Ivoire	LI	Liechtenstein	SK	Slovak Republic
CM	Cameroon	LK	Sri Lanka	SN	Senegal
CS	Czechoslovakia	LU	Luxembourg	SU	Soviet Union
CZ	Czech Republic	MC	Monaco	TD	Chad
DE	Germany	MG	Madagascar	TG	Togo
DK	Denmark	ML	Mali	UA	Ukraine
ES	Spain			US	United States of America

INTERNATIONAL SEARCH REPORT

International application No.

PCT/US 92/06617

A. CLASSIFICATION OF SUBJECT MATTER

IPC5: C12Q 1/68

According to International Patent Classification (IPC) or to both national classification and IPC

B. FIELDS SEARCHED

Minimum documentation searched (classification system followed by classification symbols)

IPC5: C12Q

Documentation searched other than minimum documentation to the extent that such documents are included in the fields searched

Electronic data base consulted during the international search (name of data base and, where practicable, search terms used)

DATA BASE: MEDLINE, BIOTECHNOLOGY ABSTRACTS

C. DOCUMENTS CONSIDERED TO BE RELEVANT

Category*	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
X	Dialog Information Services, File 154, Medline, Dialog accession no. 07171566, Medline accession no. 90078566, Venkatesan MM et al: "Use of Shigella flexneri ipaC and ipaH gene sequences for the gene- ral identification of Shigella spp. and enteroinva- sive Escherichia coli", J Clin Microbiol Dec 1989, 27 (12) p 2687-91 --	1
X	Dialog Information Services, File 154, Medline, Dialog accession no. 07673589, Medline accession no. 91192589, Cleuziat P et al: "Specific detection of Escherichia coli and Shigella species using frag- ments of genes coding for beta-glucuronidase", FEMS Microbiol Lett Nov 1990, 60 (3) p 315-22 --	1

☒ Further documents are listed in the continuation of Box C. ☒ See patent family annex.

* Special categories of cited documents:	"T" later document published after the international filing date or priority date and not in conflict with the application but cited to understand the principle or theory underlying the invention
"A" document defining the general state of the art which is not considered to be of particular relevance	"X" document of particular relevance: the claimed invention cannot be considered novel or cannot be considered to involve an inventive step when the document is taken alone
"B" earlier document but published on or after the international filing date	"Y" document of particular relevance: the claimed invention cannot be considered to involve an inventive step when the document is combined with one or more other such documents, such combination being obvious to a person skilled in the art
"L" document which may throw doubts on priority claim(s) or which is cited to establish the publication date of another citation or other special reason (as specified)	"&" document member of the same patent family
"O" document referring to an oral disclosure, use, exhibition or other means	
"P" document published prior to the international filing date but later than the priority date claimed	

Date of the actual completion of the international search	Date of mailing of the international search report
2 February 1993	25.02.93

Name and mailing address of the ISA/
European Patent Office, P.B. 5818 Patentlaan 2
NL-2280 HV Rijswijk
Tel. (+31-70) 340-2040, Tx. 31 651 epo nl,
Fax: (+31-70) 340-3016

Authorized officer
MIKAEL G:SON BERGSTRAND

INTERNATIONAL SEARCH REPORT

International application No.

PCT/US 92/06617

C (Continuation). DOCUMENTS CONSIDERED TO BE RELEVANT

Category*	Citation of document, with indication, where appropriate, of the relevant passages	Relevant to claim No.
A	EP, A2, 0357306 (INTEGRATED GENETICS, INC.), 7 March 1990 (07.03.90)	1-44
	--	
A	US, A, 4816389 (SANSONETTI ET AL), 28 March 1989 (28.03.89)	1-44
	--	

Form PCT/ISA/210 (continuation of second sheet) (July 1992)

INTERNATIONAL SEARCH REPORT

International application No.

PCT/US 92/06617

Box I Observations where certain claims were found unsearchable (Continuation of item 1 of first sheet)

This international search report has not been established in respect of certain claims under Article 17(2)(a) for the following reasons:

1. ☐ Claims Nos.:
because they relate to subject matter not required to be searched by this Authority, namely:
2. ☐ Claims Nos.:
because they relate to parts of the international application that do not comply with the prescribed requirements to such an extent that no meaningful international search can be carried out, specifically:
3. ☐ Claims Nos.:
because they are dependent claims and are not drafted in accordance with the second and third sentences of Rule 6.4(a).

Box II Observations where unity of invention is lacking (Continuation of item 2 of first sheet)

This International Searching Authority found multiple inventions in this international application, as follows:

The claimed invention relates to DNA fragments derived from the genus *Shigella* that are useful when detecting *Shigella* and/or Enteroinvasive *E. coli*. Both the problem (Detecting *Shigella* bacteria) and the solution (Detecting *Shigella* bacteria by using DNA probes derived from the *Shigella* chromosome) are known in the art. This leads to the following regrouping:

1. ☐ As all required additional search fees were timely paid by the applicant, this international search report covers all searchable claims.
2. ☒ As all searchable claims could be searched without effort justifying an additional fee, this Authority did not invite payment of any additional fee.
3. ☐ As only some of the required additional search fees were timely paid by the applicant, this international search report covers only those claims for which fees were paid, specifically claims Nos.:
4. ☐ No required additional search fees were timely paid by the applicant. Consequently, this international search report is restricted to the invention first mentioned in the claims; it is covered by claims Nos.:

Remark on Protest

- ☐ The additional search fees were accompanied by the applicant's protest.
- ☐ No protest accompanied the payment of additional search fees.

FURTHER INFORMATION CONTINUED FROM PCT/ISA/

1. Claims 1,4,7-10,17-21, and 39-44: Probes derived from the chromosome of *Shigella flexneri* and probe compositions comprising these probes.
2. Claims 2-3,5-6,12-38: Probes derived from the chromosome of *Shigella sonnei* and probe compositions comprising these probes.
3. Claims 11 and 17-26: Probes derived from the chromosome of *Shigella dysenteriae* and probe compositions comprising these probes.

INTERNATIONAL SEARCH REPORT
 Information on patent family members

08/01/93

International application No.

PCT/US 92/06617

Patent document cited in search report	Publication date	Patent family member(s)	Publication date
EP-A2- 0357306	07/03/90	JP-A- 2238899 US-A- 5084565	21/09/90 28/01/92
US-A- 4816389	28/03/89	EP-A,B- 0170584 FR-A,B- 2567541 JP-A- 61044000 US-A- 4992364	05/02/86 17/01/86 03/03/86 12/02/91